

DOVER™ | EXAM

PORTABLE LED HEADLIGHT SYSTEM



SEE BETTER
SEE FARTHER
SEE WHY


SINCE 1971
ENHANCING
YOUR VISION

Powered by
Long Island
Technology
Group

DOVER™ EXAM

HIGH-INTENSITY, PORTABLE LED EXAM HEADLIGHT WITH FIXED SPOT, POWERED BY A 3-CELL BATTERY

- Illumination Range: ~30,000 LUX at 16" (40 cm) to ~105,000 LUX at 10" (25 cm)
- Fixed Spot Diameter: 3" (7.6 cm) at a working distance of 16" (40 cm)
- Color Rendering Index: 78 CRI
- Color Temperature: 6,500 K
- Intensity Settings: 5 Levels
- 3-Cell Battery
- Battery Duration at Full Intensity: 10 hours
- Battery Expected Life Span: > 300 Charges
- LED Expected Life Span: > 10,000 hours
- Optic Weight: 0.45 oz (12.7 g)
- Battery Weight: 6 oz (170 g)



DOVER PORTABLE LED HEADLIGHT SYSTEM INCLUDES:

- Dover High-Intensity Exam Headlight
- Choice of: Fiji or Key West Headband, or Optic Clip
- Choice of One or Two 3-Cell Battery
- Charging Cradle, Power Supply
- Holster, Linkage Kit, and Carrying Case

DOVER PORTABLE LED HEADLIGHT SYSTEM PRODUCT CODES:

BFW 9030	Dover Headlight System with Fiji Headband with One 3-Cell Battery
BFW 9030B	Dover Headlight System with Fiji Headband and Two 3-Cell Batteries
BFW 9032	Dover Headlight System with Optic Clip with One 3-Cell Battery
BFW 9032B	Dover Headlight System with Optic Clip and Two 3-Cell Batteries
BFW 9034	Dover Headlight System with Key West Headband with One 3-Cell Battery
BFW 9034B	Dover Headlight System with Key West and Two 3-Cell Batteries
BFW 9035	Dover 3-Cell Battery
BFW 9036	Dover Holster
BFW 9022ACC	Dover Charger
BFW 1552C	Dover Optic Clip

EXCELLENT FOR USE IN THE FOLLOWING DEPARTMENTS

HOSPITAL	ENT GENERAL SURGERY	SURGERY CENTER	PLASTIC SURGERY
OFFICE	DENTAL EXAM ORAL/MAXILLOFACIAL		



BFWinc.com

800-717-4673 | +1-502-899-1808

445 BAXTER AVENUE | SUITE 175
LOUISVILLE, KY 40204 USA



DOVER TESTED AND CERTIFIED TO:
IEC 60601-1 Edition 3, IEC 60601-2 Edition 4
ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10),
CAN/CSA-C22.2 NO 60601-1- (2008), CE

ENGINEERED AND ASSEMBLED IN THE U.S.A. © 2020 BFW, Inc.
ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE