

Systeme de micro-caméra Full HD 3MOS

GP-UH332

Systeme de micro-caméra à haute résolution jusqu'à 1 000 lignes TV



KEY FEATURES

- Tête de micro-caméra Full HD avec capteur d'images 1/3"
- Normes sur les équipements médicaux électriques IEC60601
- Excellente reproduction des couleurs
- Sorties vidéo HDMI et HD-SDI
- Facilité de mise en œuvre dans de nombreuses applications

- Tête de micro-caméra Full HD de 1/3 de pouce
- Conforme aux normes sur les équipements médicaux électriques IEC60601*
- Excellente reproduction des couleurs
- Sorties vidéo HDMI et HD-SDI
- Facilité de mise en œuvre dans de nombreuses applications

SPECIFICATIONS

| | |
|---------------------------------|--|
| Operating Temperature | 0°C to 40°C (32°F to 104°F) |
| Power Supply | 12 V DC |
| Power Consumption | Approx. 40 W |
| Operating Humidity | 30° to 90° (without condensation) |
| Optical System | Microprism optical system |
| Horizontal Resolution | 1000 TV line standard |
| S/N | 54 dB (Gain off) |
| Standard Illumination | 2000 lx (3200K, F11) |
| Finish | Chemical-resistant coating (sail white) |
| Cable Length | 3m, 6m until the end of 2016, afterwards 6.6m, 15m |
| External Control | ø3.5 stereo mini jack |
| Image Sensor | 1/3" progressive color 3MOS |
| Video HDMI Output | HDMI x 2ch (Full HD) |
| Video SDI Output | 3G/HD-SDI x 5ch (Full HD) |
| Video Format HDMI Output | 1920 x 1080/59.94p, 1920 x 1080/50/94i, 1920 x 1080/50p, 1920 x 1080/50i |
| Video Format SDI Output | 1920 x 1080/59.94p, 1920 x 1080/50.94i, 1920 x 1080/50p, 1920 x 1080/50i |

URL: <https://business.panasonic.fr/camera-professionnelle/systeme-de-micro-camara-haute-resolution/gp-uh332>

CONTACT

Web: <https://business.panasonic.fr/camera-professionnelle/contact-us>