

ATON™ | SURGICAL

FIBER OPTIC HEADLIGHT



SEE BETTER  
SEE FARTHER  
SEE WHY

SINCE 1971  
ENHANCING  
YOUR VISION

Powered by  
Long Island  
**Technology**  
Group



**AtoN™ SURGICAL**

**HIGH-INTENSITY SURGICAL FIBER OPTIC HEADLIGHT WITH ADJUSTABLE SPOT, POWERED BY HATTERAS LED LIGHT SOURCE**

- Adjustable Iris Spot Diameter: .5"-4" (1.27 - 10.16 cm) at a working distance of 16" (40 cm)
- Color Rendering Index with Hatteras: 70 CRI
- Color Temperature with Hatteras: 5,400 K
- Standard Accompaniment: 10' Cobra Cable
- Hygienic Optic: Easy to Clean Twist-Off Feature
- Compatibility: ACMI Port
- Optic Weight: 1.5 oz (43.4g)

**HATTERAS/AtoN SYSTEM**

- Hatteras LED High-Intensity Surgical Light Source
- AtoN High-Intensity Headlight
- Choice of: Fiji or Key West Headband
- 10' Cobra Coil Cable
- Mobile Stand

**AtoN HEADLIGHT SYSTEM PRODUCT CODES:**

BFW 1200	AtoN Headlight with Fiji Headband and 10' Cobra Coil Fiber Optic Cable
BFW 1252	AtoN Headlight with Key West Headband and 10' Cobra Coil Fiber Optic Cable
BFW 1255	AtoN Headlight with Fiji Headband and 10' Standard Fiber Optic Cable
BFW 1257	AtoN Headlight with Key West Headband and 10' Standard Fiber Optic Cable
BFW 1260	AtoN Headlight with Fiji Headband and 7' Standard Fiber Optic Cable
BFW 1262	AtoN Headlight with Key West Headband and 7' Standard Fiber Optic Cable
BFW 9860	7' AtoN Standard Fiber Optic Cable
BFW 9865	10' AtoN Standard Fiber Optic Cable
BFW 9866	10' AtoN Cobra Coil Fiber Optic Cable



**EXCELLENT FOR USE IN THE FOLLOWING DEPARTMENTS**

**HOSPITAL**

- CARDIOSURGERY
- NEUROSURGERY
- ORTHOPEDIC SURGERY
- PEDIATRIC SURGERY
- SPINE SURGERY
- THORACIC SURGERY
- TRANSPLANT SURGERY
- VASCULAR SURGERY

**SURGERY CENTER**

- ORTHOPEDIC SURGERY



**BFWinc.com**

**800-717-4673 | +1-502-899-1808**

445 BAXTER AVENUE | SUITE 175

LOUISVILLE, KY 40204 USA

**ENGINEERED AND ASSEMBLED IN THE U.S.A.**

**HATTERAS TESTED AND CERTIFIED TO:** IEC 60601-1-2 Edition 4, IEC 60601-2-18 Edition 3, ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10), CAN/CSA-C22.2 NO 60601-1- (2008), CE

ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

© 2020 BFW, Inc.



Powered by  
**Long Island Technology Group**