

sensors

batteries

equipment

# Axcom: High Performance for medicine.

certified quality, 100% compliant and with up to 50% price advantage

# When it matters: **Axcom.**

The demands placed on medical equipment in hospitals and clinics are constantly increasing. Equipment must be reliable all around the clock. As a leading manufacturer of rechargeable batteries and supplier of medical accessories, Axcom makes an important contribution to the effective care of patients worldwide.

#### Certified quality without any compromises.

Reliability, endurance and longevity are Axcom's top priorities. This is why development, production and final inspection all take place in our own workshop, with quality double-checked in both production and shipping. Even the manufacturers' original batteries as well as the strictly selected accessories are carefully checked before they are delivered. All products are certified and comply 100% with the German medical products act, MPG.



- Top quality thanks to high-quality Japanese quality cells and thoroughly tested components
- Personal support by trained medical product consultants
- 100% compliance with the German medical products act, MPG
- Our own certified products and originals
- 12 month warranty for batteries in use
- Quick delivery within 24 hours
- Technical support provided by qualified team
- Immediate availability of high-demand articles through stock levels maintained by Axcom
- Significant price advantage
- Quality double-checked in both production and shipping
- Our own production, technology and quality management complying with ISO9001:2015
- Advanced technology developed in cooperation with our customers
- Tailor-made products to meet your needs

#### Axcom brand quality

Great responsibility requires top quality: Axcom's batteries are world leaders in terms of their construction, selected materials and workmanship.

We use only brand-name components to construct our batteries – from the cells to the gold-plated contacts. Our certified suppliers have been carefully selected. Every Axcom product is subject to stringent quality controls in our in-house testing laboratory. The company is also certified according to ISO 9001:2015.

#### Individual development

Axcom batteries and communication accessories are designed to your exact requirements.

Our engineers work with you from the specifications, modify standard applications and optimise products to your needs. Through our global network with leading developers and producers, Axcom products are technologically always state-ofthe-art.

We bring your product to life

#### More flexibility

Just tell us what you need - we will do it for you.

We produce small or large series, standard or special products, and supply from our automated central warehouse for just-intime to Germany or to the global offices of your company.


#### Advice and service

#### You benefit from individual project support in the team of Axcom specialists.

We advise you competently and see ourselves

as your solution partner for secure and reliable mobile energy carriers and charging technology in almost every field of application. That's why we have worked together with our customers for years on a full confidence basis.

## What Axcom offers:



## Contents

	from page		from page
ABOTT	5	HONEYWELL	22
AIR-SHIELDS	5	IMED	22
AKS / ERGOLET	5	INVIVO	23
ALARIS (IVAC)	5	JÄGER	23
AQUATEC	6	LAERDAL	23
ARCOMED	6	MAQUET	24
ARJO	6	MARC HORN	24
ATMOS	8	MARQUETTE	25
AVI/3M	8	MELA	26
BECTON DICKINSON	8	METRAX/PRIMEDIC	26
BEKA	8	MGVG/DÖRING	26
BOSCH	8	MINDRAY	26
BRAND	9	MOLIFT	27
B BRAUN	9	NELLCOR	27
BRUKER	10		28
CODAN	10		28
COROMETRIX	11	OHMEDA/BRAUN	28
CRITICARE CRITICON	11 11	PHILIPS PHYSIO CONTROL	29
CUSTOMED	11	PRIMEDIC/METRAX	30 31
DATASCOPE	12	PROTOCOL PROPAQ	31
DATASCOPE DATEX/OHMEDA	12	RADIOMETER	31
DIMEQ	13	SANITAS GE	31
DRÄGER	13	SCHELL	32
ECKLI	14	SCHILLER	32
ELECTRONICA TRENTINA	14	SCHIWA/MGVG	33
ESAOTE	14	SECA	33
FRESENIUS-VIAL/MCM	15	SIEMENS	33
GE HELLIGE/MARQUETTE	16	SIMONSEN + WEEL	34
GE MEDICAL SYSTEMS	17	SÖHNGEN	34
GRASBY	17	SPACELABS	34
GULDMAN	17	STIEGELMEYER	35
HAMILTON	18	SUNRISE MEDICAL	35
HANAU MED	18	SYNTHES	36
HANSEATISCHE INTERMED	18	TRAUTWEIN	36
HEINE	19	TRIMED	36
HEINEN + LÖWENSTEIN	19	VÖLKER	36
HELLIGE	19	VICTOR	37
HELLIGE/MARQUETTE	19	WEINMANN	37
HEWLETT PACKARD	20	WISSNER-BOSSERHOF	37
HEYMER	21	WELCH ALLYN	37
HILLROM	21	ZOLL	38
HILL-ROM/LIKO	22	CELL REPLACEMENT PATIENT LIFT	39
HORCHER	22	SPO <sub>2</sub> -SENSORS	40
HÖRMANN	22	BLOOD PRESSURE CUFF	42
		ECG-CABLES	44

ABOTT

AKS / ERGOLET

	Description / original no. Feeding pump Patrol		/pe	Voltage/capacity 6,0V / 4,5Ah	axcom product no. GB640
	Batteries compatible with			A	IR-SHIELDS
	Description / original no.	Ту	/pe	Voltage/capacity	axcom product no.
	TI 500	lea	ad	12,0V / 25,0Ah	SB51251

Batteries compatible with	•
---------------------------	---

Batteries compatible with ...

Description / original no.	Туре	Voltage/capacity	axcom product no.
patient lift Faun CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C/	lead ASING)	36,0V / 2,3Ah	MB9001
patient lift A/S Modell Maya 791 CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING CA	lead ASING)	36,0V / 2,0Ah	MBU791
lift Atlas / Atlas II / Nestor / AKS	lead	2 x12,0V / 7,2Ah	2 x MPB12-7200
lift Type Duo	lead	2 x12,0V / 4,5Ah	2 x MPB12-4500
Foldy / clino / dualo	lead	2x 12,0V/4,5Ah	2 x MPB12-4500

Batteries compatible with	ALA	RIS CAREFU	SION (IVAC)
Description / original no.	Туре	Voltage/capacity	axcom product no.
Asena CC/GS/GH/Tiva	NiMH	7,2V / 2,7Ah	MB1072
Asena GW	NiMH	7,2V / 1,3Ah	MB1078
Asena infusion pump GP	NiMH	9,6V / 2,7Ah	MB1079
syringe pump 741	lead	6,0V / 1,3Ah	GB610
syringe pump 1000-7000/Tiva	lead	6,0V / 3,4Ah	GB626

Batteries compatible with			AQUATEC
Description / original no.	Туре	Voltage/capacity	axcom product no.
patient lift Aquatec Fortuna / Beluga	lead	12,0V / 3,4Ah	MPB12-3400
Batteries compatible with			ARCOMED
Description / original no.	Туре	Voltage/capacity	axcom product no.
infusion pump Volumed VP5000 / VP5005 / VP7000	NiMH	9,6V / 2,15Ah	MB2100
syringe pump Syramed SP 6000	NiMH	9,6V / 2,15Ah	MB2100
Batteries compatible with			ARJO



Description / original no.	Туре	Voltage/capacity	axcom product no.
carendo chair/Maxi Twin lift battery type NAE0100-083 CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING CA	lead SING)	24,0V / 2,5Ah	MBU083
combilifter	lead	12,0V / 7,2Ah	GB1270A
Contoura	lead	4x 12,0V / 7,2Ah	4x GB1270
shower trolley concerto	lead	24,0V / 1,2Ah	2 x FIB12-1200

# Accessories for Arjo patient lift.

recently. Axcom offers alternatives of exceptional technology are available. Please contact our medical quality and at a considerably lower price.

Arjo products were only available as originals until Various different kinds of manual switch and charging product consultants for more information.



Product no.: SM-AR0101 Charging devices for batteries Arjo Huntleigh Tempo/Opera/Encore/ Sara Plus lift

Also available for other models. Please contact our medical product consultants.



#### Product no.: HSA-4T

Manual switch for Arjo Huntleigh Sara3000/ Sara Plus lift

Also available for other models. Please contact our medical product consultants.

	Batteries compatible with			ARJO
	Description / original no.	Туре	Voltage/capacity	axcom product no.
eplacement	patient lift Maximove ES412D/Alenti/Encom/Miranti Ma NDA0100-03 – (please send in housing) CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING		24,0V / 4,0Ah	MBU1301
	lift Maximove	lead	24V / 4Ah	MB1301
	patient lift Maximove/Ambulift/Maxilift	lead	12,0V / 4,0Ah	MPB12-4D
r al	patient lift Maximove ES412D/Alenti/Encom/Miranti Marisa type NDA0100-03 - battery insert -	lead	12,0V / 4,0Ah	DB12-4 (2x)
placement.	patient lift Alenti/Encom/Miranti/Marisa/ Sara 3000 Maximove NDA0100-03/NDA0100-20 CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING	lead GCASING)	24,0V / 4,0Ah	MBU0100
	patient lift Sara3000/Maximove/Marisa/Miranti/Bolero/ CalypsoNDA0100-03/NDA0100-20 /110212	lead	24,0V / 4,0Ah	MB0100
	patient lift Tempo/Opera/ Encore/Sara PLUS KPA0100 /110325	lead	24,0V / 6,0Ah	MB0101
	patient lift lift Tempo/Opera type KPA0100-0510	lead	24,0V / 5,0Ah	MBU0101
	patient lift Triexie/Sarita KTA 0102	lead	24,0V / 2,0Ah	MBU0102
- Martin	lift type Minstrel HMX00310.0a	lead	2x12,0V / 5,0Ah	2x YB12-500
	wall charger for Arjo Maxi Move NDA8200/HMX987		24,0V charger	SM-AR0100
-	wall charger for Arjo Tempo/ Opera/ Encore/ Sara Plus KPA8200	SPL3024,	24,0V/ 1A	SM-AR0101
ļ	battery charger for Arjo Tempo / Opera / Encore / Sara Sarita / TrixieKPA2000/KPA8200/SPL3024	Plus /	24,0V/ 1A	SM-AR0101NT
	wall charger for Arjo rechargeable battery SARA LIft 20 KTA0101-EU	00	24,0V/ 1A	SM-AR0102
	wall charger for Arjo Ambulift / Maxi Move (battery ch Maxilift	arger) /	24,0V/ 1A	SM-AR1301
HARD S. A.: Y	battery charger for Arjo rechargeable battery Ambulift Move (battery charger) / Maxilift 24V/ 6Ah	/ Maxi	24,0V/ 1A	SM-AR1301NT
)				

Batteries compatible with			ATMOS
Description / original no.	Туре	Voltage/capacity	axcom product no.
Atmolit suction pump	NiCd	12,0V / 1,3Ah	MB1096



Description / original no.	Туре	Voltage/capacity	axcom product no.
infusions pump100/110/200/200A/270/285/400A	lead	12,0V / 2,3Ah	YB12-2300
450/480/840/880/Micro 285			

#### **BECTON DICKINSON**

AVI/3M

F	

Description / original no.	Туре	Voltage/capacity	axcom product no.
Program 1	lead	6,0V / 1,2Ah	SB5611
Program 2	lead	6,0V / 3,5Ah	SB5631

	Batteries compatible with			BEKA
	Description / original no.	Туре	Voltage/capacity	axcom product no.
R.	Carlo classic shower lift	lead	12,0V / 7,2Ah	MPB12-7200

	Batteries compatible with			BOSCH
	Description / original no.	Туре	Voltage/capacity	axcom product no.
A Real Post	ECG 501 N/NV	lead	12,0V / 0,8Ah	MB1041B
- STAR	ECG 503	lead	12,0V / 1,3Ah	GB1212

В

**B**|**B**RAUN

	-	0		
		-		
-	-			
-		ALL DOOR DOOR DOOR DOOR DOOR DOOR DOOR DO	R	

Batteries compatible with			BRAND
Description / original no.	Туре	Voltage/capacity	axcom product no.
Pipetting robot	NiCd	2,4V / 1,0Ah	MB9822

Batteries compatible with ...



Description / original no.	Туре	Voltage/capacity	axcom product no.
Top-sellers made in-house: With quality cel	ls byPan	asonic:	
Infusomat fmS/P	NiCd	7,2V / 1,2Ah	MB1125P
Infusomat/ Perfusor FM	NiCd	7.2V / 2.0Ah	MB1004P
Perfusor FM (MFC) 34508210	NiCd	7,2V / 1,2Ah	MB1008P
Perfusor compact	NiCd	3,6V / 0,6Ah	MB1120
Perfusor FT/Perfusor VII - slim design -	NiCd	7,2V / 1,3Ah	MB1000
BRAUN Perfusor FT91/Perfusor F/Secura FT (triangular shape)	NiCd	7,2V / 1,2Ah	MB1002
Perfusor/Infusomat Space - original (with PIN)	NiMH	4,8V / 2,1Ah	MB1126MP-O
Perfusor/Infusomat Space CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C/	NiMH ASING)	4,8V / 2,1Ah	MBU1126
Puls Oxymeter 3800/3900	lead	8V / 3,2Ah	MB3800



Prod. no. MBU1126 NiMH battery (cell replacement / repair order) Insert for Braun Perfusor / Infusomat Space 4.8V/2.5 Ah





#### MORE CHARGING AND DISCHARGING CYCLES:

Axcom battery cells are measured individually to ensure performance and durability. Only cells with low charging resistance are used.

Batteries compatible with ...

Description / original no.	Туре	Voltage/capacity	axcom product no.
defibrillator Defigard 1002/2000/2002/semi-automatic machine Fred	NiCd	12,0V / 2,4Ah	MB2001
defibrillator MR	NiCd	12,0V / 1,3Ah	MB1029
Fred easy original lithium battery	Li	12,0V / 2,8Ah	MB2013-O
Fred original lithium battery	Li-lon	18,0V / 4,75Ah	MB2003-O
Schiller - Defigard 3002	NiCd	14,4V / 1,9Ah	MB2004



Batteries compatible with			CODAN
Description / original no.	Туре	Voltage/capacity	axcom product no.
syringe pump Argus 600	lead	6,0V / 4,5Ah	GB640
infusions pump707	NiMH	12,0V / 1,8Ah	MB407

#### Quality in the details

The longevity and quality of Axcom batteries are the result of meticulousness and perfection. The highest degree of precision, which shows in the details. That includes our durable, ultra-precise plug-in connectors



	Batteries compatible with		C	OROMETRIX	B
E	Description / original no.	Туре	Voltage/capacity	axcom product no.	С
	Medical monitor 511	lead	6,0V / 4,5Ah	GB640	
-14	Medical monitor 555	lead	12,0V / 2,2Ah	GB1219	

Batteries compatible with		
Description / original no	Type	Voltage/capa



-

Description / original no.	Туре	Voltage/capacity	axcom product no.
pulse oximeter 503	lead	6,0V / 1,3Ah	GB610

	Batteries compatible with			CRITICON
	Description / original no.	Туре	Voltage/capacity	axcom product no.
in the second se	Dinamap Compact/S/Plus	lead	12,0V / 2,0Ah	GB2312M
	Dinamap Pro Care 100/ Pro Care 200/ Pro Care 300 / V100	lead	6,0V / 3,3Ah	DB624
	Dinamap PRO 100/PRO 200/PRO 300/PRO 400	lead	12,0V / 2,0Ah	GB2312M
1	Dinamap 8100XL/8102//9300/9350/8120	lead	12,0V / 2,8Ah	GB1224C
A Contraction	Dinamap PRO 1000	NiMH	12,0V / 7,4Ah	MB1250
	multi-function monitor 8710/8720	lead	12,0V / 2,0Ah	GB2312M

Batteries compatible with			CUSTOMED
Description / original no.	Туре	Voltage/capacity	axcom product no.
Blood pressure meter - battery insert -	NiCd	4,8V / 0,6Ah	MB5011



	- ^ C 🦳	
JAI	ADU	OPE
<i>–</i>		· · · -



Description / original no.	Туре	Voltage/capacity	axcom product no.
transport monitor Accutor III - build year 96/97 -	lead	12,0V / 2,0Ah	FIB12-2000
Accutor Plus/Trio 0146-00-0069	Li-lon	11,1V / 4,8Ah	MB0241
transport monitor Passport Accutor Plus/Trio	lead	12,0V / 2,0Ah	GB2312M

Batteries compatible with ...

#### DATEX/OHMEDA



Description / original no.	Туре	Voltage/capacity	axcom product no.
defibrillator MR	NiCd	12,0V / 1,3Ah	MB1029
Puls Oxymeter 3800/3900	lead	8,0V / 3,2Ah	MB3800
Hand Held pulse oximeter Biox 3740	lead	12,0V / 2,2Ah	GB1219
monitor AS3 Light CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C	NiCd ASING)	13,2V / 1,8Ah	MBU1048
monitor S5, AS3	NiMH	2x12,0V / 4,0Ah	MB1049
Satellite Trans	lead	6,0V / 3,4Ah	GB626

	Batteries compatible with			DIMEQ	
	Description / original no.	Туре	Voltage/capacity	axcom product no.	С
	defibrillator 503	NiCd	18,0V / 0,6Ah	MB1023	
-	ECG 503 V	lead	12,0V / 1,3Ah	GB1212	D
- C					

	Batteries compatible with			DRÄGER
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	Babylog	lead	6,0V / 4,5Ah	GB640
	defibrillator Cardiolog 2000	NiCd	12,0V / 1,9Ah	MB2000A
	Evita 4/Evita 2 dura external	lead	12,0V / 17,0Ah	GB1217 (2x)
	Evita 4/Evita 2 dura internal	lead	12,0V / 3,5Ah	SB5131 (2x)
	Julian anesthesia-respirator	lead	12,0V / 7,2Ah	GB1270A (2x)
multipower interior network of service App 1 200	monitor Infinity SC6002XL/SC6802XL/SC7000XL/ SC9000XL	lead	12,0V / 2,0Ah	GB2012M
	Oxylog 2000	NiCd	7,2V / 1,0Ah	MB1051
	Oxylog3000 (2M86732)	NIMH	12V / 3,8Ah	MB1230
Tel	Oxylog 3000 original	Li-Ion	11,1V / 4,0Ah	MB1230-O
	monitor Dräger Infinity Gamma/Delta Gamma XL+XXL/SC6002XL/SC6802XL/DeltaXL	Li-lon	14,4V / 4,0Ah	MB1057
ALCONING ALCONING	monitor Dräger Infinity Gamma/Delta Gamma XL+XXL/SC6002XL/SC6802XL/DeltaXL	Li-Ion	14,4V / 6,8Ah	MB1057-68
Panisono Alexandro Laboratoria Alexandro Carlos and Alexandro Alexandro Alexandro Alexandro A	Oxipack 2000+C787	lead	6,0V / 1,3Ah	GB610
1	Oxysat 2	lead	6,0V / 2,8Ah	GB624



Batteries compatible with				ECKLI
Description / original no.	Ту	pe	Voltage/capacity	axcom product no.
Suction Pump EMS	lea	ad	12,0V / 7,2Ah	GB1270A



Batteries compatible with	ELECTRONICA TRENTINA			
Description / original no.		Туре	Voltage/capacity	axcom product no.
Cardiette 3 channel ECG		lead	12,0V / 0,8Ah	MPB12-800

	Batteries compatible with			ESAOTE
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	Mylab One	Li-Ion	14,4V / 6,6Ah	MB4401
-	Hamilton respirator C2 / Esaote MyLab Five BK Ultrasound Flex Focus 500	Li-Ion	14,4V / 6,6Ah	MB2024

**TOP SELLER** from **axcom** – just one example from our comprehensive product portfolio:



Ni-MH battery for

ALARIS ASENA CC/GS/GH/TIVA

7.2V/2.7Ah (original no. 1000SP01122)

- with declaration of conformity - structurally identical to original, fully compatible

Prod. no. MB1072

With axcom you always benefit from:

- tested quality and ultimate product safety engineered in our own facilities

- a comprehensive program of testing in German hospitals

- extensive service and expert advice provided by medical product consultants

axcom 14

#### < Contents

Batteries compatible with ...







Description / original no.	Туре	Voltage/capacity	axcom product no.
Feeding pump Applix	NiMH	4,8V / 1,25Ah	MB1034
infusions pump Inca/ST	lead	12,0V / 3,4Ah	GB1226
infusions pump Volumed VP 7000	NiMH	9,6V / 2,15Ah	MB2100
injectiomat C/CP	lead	6,0V / 1,3Ah	GB610
injectiomat C/CP - with printed solder tail -	NiCd	1,2V / 1,7Ah	M12TBF
syringe pump injectiomat Agilia	NiMH	6,0V / 1,9Ah	MB3639
syringe pump Injektomat Agilia - (original )	NiMH	6,0V / 1,9Ah	MB3639-O
infusions pump Volumat Agilia D	Li-lon	7,2V / 2,2Ah	MB1790
infusions pump Volumed VP2000/VP5000 -	NiCd	12,0V / 1,8Ah	MB1037-18
injectiomat C/CP - high-temperature-battery -	NiCd	1,2V / 1,2Ah	M12HTTBF
injectiomat 2000/Pilot A/A2/C	lead	6,0V / 1,3Ah	GB610
MCM 404/504/Volumed VP 2000/5000	NiCd	12,0V / 1,6Ah	MB404
MCM 440/550/OPTIMA	NiMH	6,0V / 3,0Ah	MB440

**TOP SELLER** from **axcom** – just one example from our comprehensive product portfolio:



#### With axcom you always benefit from:

- tested quality and ultimate product safety engineered in our own facilities

- a comprehensive program of testing in German hospitals

- extensive service and expert advice provided by medical product consultants

#### Li-lon battery for

#### **DRÄGER PATIENT MONITORS**

14.4V/4.2AH Dimensions: 145 x 65 x 20 mm Fits: SIEMENS SC6002XL/SC6802XL DRÄGER Infinity Gamma, Gamma XL, Gamma XXL, Delta and Delta XL (original no. MS14490)

#### - with declaration of conformity

- structurally identical to original, fully compatible

Prod. no. MB1057

F

FRESENIUS-VIAL/MCM

Batteries compatible with ...



Description / original no.	Туре	Voltage/capacity	axcom product no.
Orchestra modul DPS/MVP	lead	6,0V / 1,3Ah	GB610
Sondomat Feeding pump	NiCd	1,2V / 0,6Ah	Z-N600AAC (4x)
Programm I	lead	6V / 1,2Ah	SB5611
Programm II	lead	6,0V / 3,4Ah	GB626
syringe pump Pilot/1	lead	6,0V / 1,3Ah	GB610
syringe pump Pilot/A2	lead	6,0V / 1,3Ah	GB610

FRESENIUS-VIAL/MCM

GE HELLIGE

_				
		-		D
	10	1911	14	1
		C.C.		
1	-	-		4









Description / original no.	Туре	Voltage/capacity	axcom product no.
Cardiotest EK 53/56	NiCd	21,6V / 1,3Ah	MB1027
defibrillator Responder 3000 original	NiCd	12,0V / 2,0Ah	MB2009-O
defibrillator SCP 844 - Akku Einsatz CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C	NiCd ASING)	12,0V / 1,3Ah	MBU2002
defibrillator SCP 910/913/915/922/Cardioserv	NiCd	12,0V / 1,2Ah	MB1030
Defiport SCP 912/SCP 840/Responder 1000/1100/2000	NiCd	18,0V / 0,6Ah	MB1023
Defiport/Defiscope M	NiCd	12,0V / 1,3Ah	MB1029
Eagle monitor 1000/1006/1008/1009	NiCd	12,0V / 2,0Ah	MB2008
ECG 3500/MAC 5000/MAC5500 (original ) 9007700-001	NiMH	1,0V / 3,5Ah	MB1208
EK 33/43/403	lead	6,0V / 3,5Ah	SB5631
ECG CardioSmart/MAC1200	NiCd	18,0V / 1,6Ah	MB1019
ECG MicroSmart	lead	12,0V / 1,3Ah	GB1212
MAC1200 Backup battery	Li-Me	3V	BUMB1019

axcom 16

	Batteries compatible with			GE HELLIGE
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	MAC1600 ECG - 2035701-001, 2032095-001 (original )	Li-Ion	14,8V / 2,3Ah	MB3808
-	MAC500	lead	12,0V / 1,3Ah	GB1212
	monitor Dash 2000	NiCd	12V / 3,0 Ah	MB1206
	monitor Dash 3000/4000	Li-lon	11,1V / 4,5 Ah	MB1201
	patient monitoring solar	lead	12,0V / 7,2 Ah	GB1270
	SMP 300/320 transport monitor	lead	12,0V / 2,0 Ah	GB2312M

-	-	Descr
05:8CE	AX	
.) 10 (10.0		

Batteries compatible with ...

Description / original no.	Туре	Voltage/capacity	axcom product no.
monitor Dash 2000	NiCd	12,0V / 3,0Ah	MB1206
monitor Dash 3000/4000/5000 F-FM/monitor B30/transport pro/SM201-6 with capacity display	Li-lon	11,1V / 4,5Ah	MB1201
monitor Vital Sign Pro 400	lead	12,0V / 2,0Ah	GB2312M

Batteries compatible with			GRASBY
Description / original no.	Туре	Voltage/capacity	axcom product no.
infusions pump2000/3000er series	lead	2,0V / 2,5Ah	GL0810-4 (3x)

Batteries compatible with			GULDMANN
Description / original no.	Туре	Voltage/capacity	axcom product no.
lift Model Samson / Samson III original -No.: 933011/933001	lead	4 x 12,0V / 2,2Ah	4 x MPB12-2200



F

G

GE MEDICAL SYSTEMS



NEW

#### QUALITY THAT REDUCES COSTS

Axcom only uses safe technologies that guarantee long life-cycles, which reduces the overall cost of a battery.

Batteries compatible with			HAMILTON
Description / original no.	Тур	e Voltage/cap	pacity axcom product no.
respirator Raphael	leac	12,0V / 2,3A	Ah YB12-2300
respirator C2/Esaote MyLab Five - NL2024	Li-lo	on 14,4V / 6,6A	h MB2024

Batteries compatible with		ŀ	IANAU MED
Description / original no.	Туре	Voltage/capacity	axcom product no.
Medap secretion suction device P7020 portable suction device 0789520200	lead	12,0V / 2,0Ah	SB5121
OP stand lamp blue 30 hospital	lead	12,0V / 7,2Ah	GB1270 (2x)
OP stand lamp blue 80 hospital	lead	12,0V / 38,0Ah	GB12380VDS (2x)

Batteries compatible with	HANSEATISCHE INTERMED		
Description / original no.	Туре	Voltage/capacity	axcom product no.
Oxiplett	lead	12,0/2,0Ah	SB5121



and a second

	Batteries compatible with			HEINE
	Description / original no.	Туре	Voltage/capacity	axcom product no.
Constra Col	M2Z-X-01.89.013, X-001.99.305 original	NiMH	2,5V / 0,75Ah	MB7010
101 92 30 2.84 C	S2Z-X-01.99.333	NiMH	2,5V / 0,65Ah	MS721-0
The second second	S3Z-X-02.99.314	NiMH	3,5V / 0,65Ah	MS722-0
MOZ LINELIS BARRIES - ST	X-02.99.106, X-002.99.382 original	NiMH	3,6V / 0,77Ah	MS106-O
inome sav	X-02.99.382 (original ) (ALTE NR. X-02.99.380)	NiMH	3,5V / 1,0Ah	MB7015

HEINEN	& L	ÖWE	NST	EIN

	Batteries compatible with			HEINEN & LO	OWENSTEIN
	Description / original no.	Т	уре	Voltage/capacity	axcom product no.
T HILLING	Oxymonn II	le	ead	12,0V / 0,8Ah	MB1041B
and and a	Galileo respirator	le	ead	12,0V / 4,0Ah	MPB12-4500

	Batteries compatible with		HELLIGE/N	IARQUETTE
	Description / original no.	Туре	Voltage/capacity	axcom product no.
Ь	Cardiotest EK 53/56	NiCd	21,6V / 1,3Ah	MB1027
	defibrillator Responder 3000 original	NiCd	12,0V / 2,0Ah	MB2009-O
[	defibrillator SCP 844 - Akku Einsatz CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C	NiCd ASING)	12,0V / 1,3Ah	MBU2002
	defibrillator SCP 910/913/915/922/Cardioserv	NiCd	12,0V / 1,2Ah	MB1030
)	Defiport SCP 912/SCP 840/Responder 1000/1100/2000	NiCd	18,0V / 0,6Ah	MB1023
	Defiport/Defiscope M	NiCd	12,0V / 1,3Ah	MB1029
	Eagle monitor 1000/1006/1008/1009	NiCd	12,0V / 2,0Ah	MB2008









	Batteries compatible with		HELLIGE/M	1ARQUETTE
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	ECG 3500/MAC 5000/MAC5500 (original ) 9007700-001	NiMH	18,0V / 3,5Ah	MB1208
	EK 33/43/403	lead	6,0V / 3,5Ah	SB5631
	ECG CardioSmart/MAC1200	NiCd	18,0V / 1,6Ah	MB1019
-	ECG MicroSmart	lead	12,0V / 1,3Ah	GB1212
	MAC1200 Backup battery	Li-Me	3V	BUMB1019
T an	MAC1600 ECG - 2035701-001, 2032095-001 (original )	Li-lon	14,8V / 2,3Ah	MB3808
-	MAC500	lead	12V / 1,3Ah	GB1212
still AX	monitor Dash 2000	NiCd	12V / 3,0 Ah	MB1206
1	monitor Dash 3000/4000	NiCd	11,1V / 4,5 Ah	MB1201
-1	patient monitoring Solar	lead	12,0V / 7,2 Ah	GB1270
	SMP 300/320 transport monitor	lead	12,0V / 2,0 Ah	GB2312M

	Batteries compatible with		HEWLETT	PACKARD
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	Codemaster XL+/M1722A+B/M1723A+B/M1724A/M1758A	lead	12,0V / 4,5Ah	MB1100
	monitor/defibrillator 43100A/43110A/43120A/43130A/ 43200A	lead	12,0V / 2,5Ah	GL43100
N.S.	monitor M 3046 A	NiMH	12,0V / 3,8Ah	MBDR36MT
TY	Pagewriter 100/200/300 i - series/1771A - without connector -	lead	6,0V / 7,2Ah	GB665
	transport monitor 40488A/M 1275/12771	lead	12,0V / 2,0Ah	GB2312M

	Batteries compatible with			HEYMER
	Description / original no.	Туре	Voltage/capacity	axcom product no.
O.E.	lift Alka / Hewo	lead	2 x 12,0V / 7,2Ah	2 x MPB12-7200
	Batteries compatible with			HILL ROM
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	Avantguard 1200 96 x 49 x 60 mm	lead	24/0,8Ah	MB3803
001	Century CC/891/893 Criticare Bed Scale ups 134 x 34 x 66 mm	lead	6,0V / 3,4Ah	MB3643
	Evolution MA2B/MA8A/MA8B/CB12, 2 batteries, 97 x 43 x 53 mm	lead	24,0V / 1,2Ah	MB3644
N. N.	Evolution CS156/LI156A/I156B/LI156C/LI156D/ LE156E/ME1A/MA3A/CB18, 2 batteries, je 97 x 43 x 53 mm	lead	24,0V / 1,2Ah	MB3645
OF	Total Care, Affinity III Beds, 2 batteries, 151 x 65 x 102 mm	lead	24,0V / 7,2Ah	MB3646
1 al	Versa Care, 2 batteries, 151 x 98 x 100 mm	lead	24,0V / 12,0Ah	MB3647
	Total Care Bed Intellidrive, 2 batteries	lead	24,0V / 17,0Ah	MB3648

**TOP SELLER** from **axcom** – just one example from our comprehensive product portfolio:



With axcom you always benefit from:

- tested quality and ultimate product safety engineered in our own facilities

- a comprehensive program of testing in German hospitals

- extensive service and expert advice provided by medical product consultants

Li-lon battery for VOLUMAT AGILIA D INFUSION PUMP 7.2V/2.2Ah

(original no. 179033-R2)

- with declaration of conformity

- structurally identical to original, fully compatible

Prod. no. MB1790

Н

#### HILL-ROM/LIKO



Description / original no.	Туре	Voltage/capacity	axcom product no.
lift Golvo 7000es / 7007 / Sabina / Viking M / Sabina II / Primo ES		2x 12,0V / 2,9Ah (PB)	2x MPB12-2,9
MIR / Viking M24 LIKO / Golvo 8000/ 8008 Type 2006107 / Linak BAL20001-02 CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING		24,0V / 2,5Ah (NiMH) 24,0V / 2,5Ah (NiMH)	MBU9001 MBU8000

Batteries compatible with ...

#### HORCHER



Description / original no.	Туре	Voltage/capacity	axcom product no.
patient lift Lexa / raisa / Diana /caesar / Ekopack origir L103-1020	nal -Notead	24,0V / 2,9Ah	MBU1060
CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING	G CASING)		
patient lift Nina / Mona original -No.: N206-1012	lead	36,0V / 2,0Ah	MB9001
CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING	G CASING)		

Batteries compatible with			HÖRMANN
Description / original no.	Туре	Voltage/capacity	axcom product no.
Bioset 3600	NiCd	9,6V / 2,0Ah	MB3600

	Batteries compatible with		НО	NEYWELL
	Description / original no.	Туре	Voltage/capacity	axcom product no.
7	Invivo Millennia	lead	12,0V / 2,0Ah	GB2312M

Batteries compatible with			IMED
Description / original no.	Туре	Voltage/capacity	axcom product no.
infusions pump Gemini PC1	lead	6,0V / 7,2Ah	GB665
infusions pump Gemini PC2 800/900-er series	lead	6,0V / 1,2Ah	SB5611
infusions pump Gemini PC4	lead	12,0V / 3,4Ah	GB1226



Batteries compatible with			INVIVO
Description / original no.	Туре	Voltage/capacity	axcom product no.
compact monitor OmniTrak 3150	lead	12,0V / 7,2Ah	GB1270

154	-	-	1	
	~			
		100	0	

Batteries compatible with			JÄGER
Description / original no.	Туре	Voltage/capacity	axcom product no.
Apnoscreen - without connector -	NiCd	6,0V / 3,0Ah	MB1064



	-	

уре	Voltage/capacity	axcom product no.	1.1
iCd	2x 6,0V / 1,9Ah	STR5 (1x)	
ad	12,0V / 2,0Ah	GB2312M	L
ad	12,0V / 2,0Ah	GB2012M	
i-Mn	12,0V / 4,2Ah	MB900	
i-Mn	12,0V / 4,2Ah	MB900-O	

**I AFRDAI** 



Description / original no.	Туре	Voltage/capacity	axcom product
Suction pump (2 x 6V-NiCd Akku)	NiCd	2x 6,0V / 1,9Ah	STR5 (1x)
Suction pump LSU / LSU 4000	lead	12,0V / 2,0Ah	GB2312M
Suction pump LPSU	lead	12,0V / 2,0Ah	GB2012M
defibrillator Heartstart FR2 M3863A/M3840A	Li-Mn	12,0V / 4,2Ah	MB900
defibrillator Heartstart FR2 M3863A original	Li-Mn	12,0V / 4,2Ah	MB900-O
defibrillator HeartStart 1000	lead	12,0V / 7,2Ah	GB1270
defibrillator HeartStart 4000 type M5516A	lead	12,0V / 2,0Ah	GB2312M
defibrillator/monitor Heartstart MRx original M3538A	Li-lon	14,4V / 3,0Ah	MB3538A-O
defibrillator/monitor Heartstart MRx NEW! ORIGINAL AXCOM	Li-Ion	14,4V /6,3Ah	MB3538A
PHILIPS/Heart S tart FR3 - (original )	Li-lon	12V / 4,7Ah	MB3002-O

Annual Control of Cont
Ministration
Anticipation of the second sec
Contraction of the second seco
and the second

Batteries compatible with			MAQUET
Description / original no.	Туре	Voltage/capacity	axcom product no.
Medap P4020/P4030/P7020/P7030 Suction Pump	lead	12,0V / 2,0Ah	SB5121
Maquet Servo I	NiMH	12,0V / 3,8Ah	MB1127
Maquet Servo I original	NiMH	12,0V / 4,0Ah	MB1127-O
OP column 115001A1	lead	12,0V / 8,0Ah	MB2020 (4x)
OP column 1150.01.C2	lead	12,0V / 14,0Ah	SBC51152 (4x)
OP column 1150.02CO	lead	12,0V / 14,0Ah	SBC51152 (4x)
OP table	lead	12,0V / 12,0Ah	GB12120A
OP table 1211.02CO	lead	12,0V / 10,0Ah	SB51101 (4x)
OP table 1140.01D	lead	12,0V / 6,5Ah	SB5161 (4x)
OP table remote control Stierlen	NiCd	10,8V / 0,6Ah	MS9180

Batteries compatible with	MARC HOR		
Description / original no.	Туре	Voltage/capacity	axcom product no.
stainless steel lift type ESMH01001-00409	lead	2 x 12,0V / 5,0Ah	2 x MPB12-5C



#### **ALWAYS THERE FOR YOU**

Question or suggestions? Looking for a product we don't have? At Axcom, service and advice are our top priority. We look forward to hearing from you: +49 (0) 21 54 - 48 38-0 / info@axcom.de

	Batteries compatible with			MARQUETTE
	Description / original no.	Туре	Voltage/capacity	axcom product no.
11ho	Cardiotest EK 53/56	NiCd	21,6V / 1,3Ah	MB1027
	defibrillator Responder 3000 original	NiCd	12,0V / 2,0Ah	MB2009-O
	defibrillator SCP 844 - rechargeable battery insert CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C	NiCd ASING)	12,0V / 1,3Ah	MBU2002
and a second sec	defibrillator SCP 910/913/915/922/Cardioserv	NiCd	12,0V / 1,2Ah	MB1030
(TOP)	Defiport SCP 912/SCP 840/Responder 1000/1100/2000	NiCd	18,0V / 0,6Ah	MB1023
	Defiport/Defiscope M	NiCd	12,0V / 1,3Ah	MB1029
	Eagle monitor 1000/1006/1008/1009	NiCd	12,0V / 2,0Ah	MB2008
	ECG 3500/MAC 5000/MAC5500 (original ) 9007700-001	NiMH	8,0V / 3,5Ah	MB1208
	EK 33/43/403	lead	6,0V / 3,5Ah	SB5631
	ECG CardioSmart/MAC1200	NiCd	18,0V / 1,6Ah	MB1019
-	ECG MicroSmart	lead	12,0V / 1,3Ah	GB1212
	MAC1200 Backup battery	Li-Me	2,5V / 0,75Ah	BUMB1019
JE	MAC1600 ECG - 2035701-001, 2032095-001 (original )	Li-Ion	14,8V / 2,3Ah	MB3808
- ORN	MAC500	lead	18,0V / 1,6Ah	GB1212
Sherry AX	monitor Dash 2000	NiCd	12V / 3,0 Ah	MB1206
P	monitor Dash 3000/4000	Li-Ion	11,1V / 4,5 Ah	MB1201
	patient monitoring solar	lead	12,0V / 7,2 Ah	GB1270
	SMP 300/320 transport monitor	lead	12,0V / 2,0 Ah	GB2312M

Batteries compatible with			MELA
Description / original no.	Туре	Voltage/capacity	axcom product no.
Vitacard	lead	12,0V / 2,2Ah	GB1219



#### METRAX/PRIMEDIC



Description / original no.	Туре	Voltage/capacity	axcom product no.
defibrillator B M110	NiCd	14,4V / 1,5Ah	MB1110
Defi HeartSave 6/HS6/HS6S/M250-97641	Li-Ion	15,0V / 7,2Ah	MB250-O
defibrillator M240 CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C	NiCd <mark>ASING</mark> )	14,4V / 1,7Ah	MBU240

#### Batteries compatible with ...

# Sour A

Batteries compatible with		MGVC	G/DÖRING
Description / original no.	Туре	Voltage/capacity	axcom product no.
Combimat IP 83/85	NiCd	6,0V / 3,0Ah	MB1011
Combimat IP 85-1+2/IP 83-2/Combimat 2000	NiCd	7,2V / 1,2Ah	MB1065
Combimat IP 83/85/Schiwamatic 7000/9000 - 6V/3AH	NiCd	6,0V / 3,0Ah	MB1011
Combimat 2000/IP85 (1+2)/CS01	NiCd	7,2V / 1,2Ah	MB1065
Combimat CS 03 / CS 06 -(large version)	NiCd	9,6V / 1,2Ah	MB1065L

	Description / origina
	VS800 monitor
	Mec 1200 monitor
NEW	IMEC12, IMEC10, IME VS600, VS90

Batteries compatible with			MINDRAY
Description / original no.	Туре	Voltage/capacity	axcom product no.
VS800 monitor	lead	12,0V / 2,0Ah	GB2312M
Mec 1200 monitor	lead	12,0V / 2,0Ah	GB2312M
IMEC12, IMEC10, IMEC8, IPM10, IPM12, IPM8, VS600, VS90	Li-Ion	11,1V / 4,5Ah	MB6008

Batteries compatible with			MOLIFT
Description / original no.	Туре	Voltage/capacity	axcom product no.
ceiling lift Type Hi-Trac	lead	2 x 12,0V / 2,2Ah (PB)	2x MPB12-2200
Molift Power Pac type 1341000	NiMH	26,4V / 3,0Ah	MB9001
Molift Power Pac type 0540000	NiCd	14,4V / 1,9Ah	MB9001
Molift Power Pac type 0541000-2	NiMH	14,4V / 3,0Ah	MB9001
Molift Power Pac type 1221115	NiCd	26.4V / 3.0Ah	MB9001















MM-GS

Batteries compatible with			NELLCOR
Description / original no.	Туре	Voltage/capacity	axcom product no.
Invos cerebral oximeter	lead	12,0V / 3,5Ah	SB5131
capnometer/pulse oximeter Type NPB70/NPB75 – CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDIN	NiMH IG CASING)	7,2V / 2,5Ah	MBU2010
patient monitor E 300	lead	12,0V / 2,0Ah	GB2312M
pulse oximeter / monitoring N-180 / N-185	lead	12,0V / 2,0Ah	GB2312M
pulse oximeter / monitoring N-190 / N-290	lead	6,0V / 4,5Ah	MPB6-4500
pulse oximeter / monitoring N-195 / N-395 / N-595	lead	6,0V / 4,5Ah	MPB6-4500
pulse oximeter / monitoring N-200	lead	6,0V / 2,8Ah	GB624
pulse oximeter / monitoring N-550	lead	12,0V / 1,3Ah	GB1212
pulse oximeter / monitoring N-2000 / 2500	lead	6,0V / 4,5Ah	GB640
pulse oximeter / monitoring N-3910 / N-3920	lead	6,0V / 4,5Ah	MPB6-4500
pulse oximeter N-560 OXIMAX	NiMH	9,6V / 4,0Ah	MB2560
pulse oximeter/monitoring N-6000	lead	12,0V / 3,4Ah	GB1226
Pulsoxy Novametrix 505/515	lead	12,0V / 2,2Ah	GB1219
pulse oximeter NPB4000	lead	6,0V / 8,0Ah	GL0859-12
	lead	6,0V / 4,5Ah	FIB6-4500

Ν

	Batteries compatible with		NIHC	N KOHDEN
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	defibrillator Cardiolife TEC5521/TEC5531/TEC7221/TEC7621 TEC7631/TEC7721/TEC7731 – original type X065/NKB-301V	NiMH	12,0V / 2,8Ah	MB2007-O
	defibrillator Cardiolife TEC5521/TEC5531/TEC7221/TEC7621 TEC7631/TEC7721/TEC7731	NiMH	12,0V / 2,8Ah	MB2007
	defibrillator ECG8823G/9010/9020/9020K/9022/9110 9130/9132 Type X041A	lead	12,0V / 2,0Ah	GB2312M
	defibrillator Lifescope 6 / OEC 6102A / 6105A	lead	12,0V / 2,0Ah	MB2012
	defibrillator 8340	lead	6,0V / 12,0Ah	GB6200 (2x)
Constant of the second	defibrillator SB901D	NiMH	12,0V / 1,95Ah	MB9010-O
	ECG84xx/ECG882x/ECG8830 type X041A	lead	12,0V / 2,0Ah	GB2312M
1 == )	monitor BSM2300 Lifescope i/BSM2301 type X062/YS-076P5	NiMH	12,0V / 3,7Ah	MB2005
-	monitor BSM2300 Lifescope i/BSM2301 type X062/YS-076P5	NiMH	9,6V / 3,7Ah	MB2015-O

Batteries compatible with		NONIN	MEDICAL
Description / original no.	Туре	Voltage/capacity	axcom product no.
pulse oximeter 8600/8604/8700/8800	NiCd	6,0V / 1,3Ah	MS8400

Batteries compatible with		OHMED	A/BRAUN
Description / original no.	Туре	Voltage/capacity	axcom product no.
pulse oximeter 370/Biox 3740	lead	12,0V / 2,2Ah	GB1219
Datex Puls Oxymeter 3800/3900	lead	8V / 3,2Ah	MB3800

Sea	

defibrillator Heartstart HS1/F M5070A - (original)

Description / original no.
















Ν

0

Ρ

PHILIPS

axcom product no.

MB5070A-O

Voltage/capacity

9.0V / 4.2Ah

Li-Mn

Lifepak 5/10 defibrillator/250 monitor

Lifepak 12 defibrillator – lead-original rechargeable

Description / original no.

Lifepak CR+ – original

Lifepak 6/7 defibrillator

- with Hawker cells -

Lifepak 9 defibrillator

battery 11141-000028

Lifepak 9 defibrillator LP9H

320-2176-001



axcom product no.

LPCRP-O

LP5FP-13

LP6200

LP9H

LP9

LP12-0

Voltage/capacity

12,0V / 1,5Ah

12,0V / 1,3Ah

12,0V / 1,8Ah

16,0V / 2,5Ah

16,0V / 2,5Ah

12,0V / 2,5Ah

Туре

Li-Me

NiCd

NiCd

lead


lead

lead











#### **TESTED FOR 100% PERFORMANCE**

Axcom batteries are tested in close collaboration with public and teaching hospitals. They comply with the stringent demands of the German Medical Products Act (MPG).



#### A RELIABLE PARTNER TO THOSE WHO ARE NEEDED

Question or suggestions? Looking for a product we don't have? At Axcom, service and advice are our top priority. We look forward to hearing from you: +49 (0) 21 54 - 48 38-0 / info@axcom.de

Batteries compatible with ...

#### PRIMEDIC/METRAX



Description / original no.	Туре	Voltage/capacity	axcom product no.
defibrillator B M110	NiCd	14,4V / 1,5Ah	MB1110
Defi HeartSave 6/HS6/HS6S/M250-97641	Li-Me	15,0V / 7,2Ah	MB250-O
defibrillator M240 CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING C	NiCd ASING)	14,4V / 1,7Ah	MBU240

Batteries compatible with ...

#### PROTOCOL PROPAQ



Description / original no.	Туре	Voltage/capacity	axcom product no.
Propaq monitor 102 / 104 / 106 / 204 / 206 - with SP02 display -	lead	8,0V / 5,0Ah	MB6082

	Batteries compatible with		R	ADIOMETER	
-	Description / original no.	Туре	Voltage/capacity	axcom product no.	
Panasonic Materials Materi	pulse oximeter OXY 3	lead	6,0V / 1,3Ah	GB610 (2x)	
	TCM 3/Tina	lead	12,0V / 3,4Ah	GB1226	Ρ
-	device TCM	lead	12,0V / 2,0Ah	GB2312M	
	TCM TOSCA monitor	lead	12,0V/2,3Ah	YB12-3200	R

S

Batteries compatible with			SANITAS GE
Description / original no.	Туре	Voltage/capacity	axcom product no.
Healthcare Pro Care B50/B40/B30/B20 monitor 2044978-004 (original )	Li-Ion	11,1V / 3,5Ah	MB1250-O

Batteries compatible with ...

AT1 ECG/spiro meter Spirovit SP1 AT101

spiro meter SP200/AT102 Cardiovit

AT2/AT2Plus/AT6/AT60/CM8/Argus-monitor ECG CS100

Description / original no.

AT3 ECG

#### SCHELL INDUSTRIES

Voltage/capacity

12,0V / 1,2Ah

12,0V / 2,2Ah

9,6V / 0,7Ah



Description / original no.	Туре	Voltage/capacity	axcom product no.
hospital bed Mod.D5912.1	lead	24,0V / 1,2Ah	MB3650
patient bed Modell J7433.1, 2 Stück	lead	24,0V / 1,2Ah	MB3649

Туре

lead

lead

NiCd

Li-lon

7,4V / 4,2Ah

A CONTRACTOR OF
FDK
FDK











AT102 plus

SCHILLER

YB12-1200

GB1219

MB1043

CR14250P

CR14250P

MB2001

MB2003-O

MB2004

MB1010LI-O

MB2013-O

MB2104

GB1226

MB1102

Batteries compatible with		SCHIV	VA/MGVG



Description / original no.	Туре	Voltage/capacity	axcom product no.
Schiwamatic 7000/9000	NiCd	6,0V / 3,0Ah	MB1011



Description / original no.	Туре	Voltage/capacity	axcom product no.
infant scales 717	NiMH	7,2V / 1,2Ah	MB1091
weighing scales 771/772/909/921/944/959	NiMH	7,2V / 0,7Ah	MB1099

#### Batteries compatible with ...

# 

Description / original no.	Туре	Voltage/capacity	axcom product no.
Cardiostat 31/31S	lead	12,0V / 2,2Ah	GB1219
monitor SC 5000/6000/7000	lead	12,0V / 2,0Ah	GB2012M
monitor SC 6002XL/6802XL converted version of SC5000/6000	Li-lon	14,4V / 4,0Ah	MB1057
SC6002XL/SC6802XL/Dräger Infinity Gamma/Delta MS30502/MS31385	Li-lon	14,4V / 6,8Ah	MB1057-68
monitor SC 7000 XL/9000 XL external	lead	12,0V / 2,0Ah	GB2312M
monitor SC 7000 XL/9000 XL internal	lead	12,0V / 3,4Ah	GB1226
monitor SC 8000	lead	12,0V / 2,2Ah	GB1219
monitor SC 9000/7000	NiCd	2x6,0V / 1,4Ah	MB1059

SECA

SIEMENS



SIEMENS
e Voltage/capacity axcom product no.
6,0V / 3,5Ah <b>SB5631 (2x)</b>
6,0V / 4,5Ah <b>GB640</b>
12,0V / 2,0Ah <b>SB5121</b>
12,0V / 2,2Ah <b>GB1219</b>



Batteries compatible with	SIMONSEN & WEE		EN & WEEL
Description / original no.	Туре	Voltage/capacity	axcom product no.
defibrillator Cardio Aid	NiCd	12,0V / 1,9Ah	MB2000A
defibrillator DMS700/DMS730/DMS750	NiCd	12,0V / 3,0Ah	MB1007
Pulsoxy 503	lead	6,0V / 1,3Ah	GB610

#### SÖHNGEN



Batteries compatible with				SÖHNGEN
Description / original no.	Ту	pe	Voltage/capacity	axcom product no.
Suction Pump (without side clips)	lea	ad	12,0V / 3,4Ah	GB1226

	Batteries compatible with			SPACELABS
-5 1	Description / original no.	Туре	Voltage/capacity	axcom product no.
EST.	monitor 1050	NiMH	12,0V / 2,4Ah	MB2051
multismus mit	PC express monitor 90305/90306/90308	lead	12,0V / 2,2Ah	GB2312M



#### **GUARANTEED RELIABILITY**

Batteries compatible with ...

Axcom backs up the quality of its batteries with a wide range of warranty services.

and the second s	
state ce	

Description / original no. Туре Voltage/capacity Vivendo EM1/EM3 lead 24,0V / 1,2Ah MB3651 2 batteries, 97 x 43 x 53 mm Seta, Type Linak BA20003-00 lead 24,0V / 1,2Ah MBU3804S Splash-proof seal, Epostyle CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING CASING) BA1812-0500-000/BA1811-24,0V / 1,2Ah MBU3804 lead 0900-000 BA18013-0, Epostyle CELL REPLACEMENT (PLEASE SEND IN BATTERY INCLUDING CASING) 24,0V / 1,2Ah MB3652 Avanta, Avanta Plus, Avanta lead Classic, EM4 2 batteries, 97 x 48,5 x 55 mm Vitano/ AG300/ Sicuro Plus lead 24,0V / 1,2Ah MB3653 4 batteries, 97 x 25 x 57 mm



Batter	ies compatible with		SUNRISE	MEDICAL
Descr	iption / original no.	Туре	Voltage/capacity	axcom product no.
patien	t lift Oxford Presence	lead	2 x 12,0V / 2,9Ah	2 x MPB12-2900
patien	t lift Oxford Classic 150 E	lead lead lead	2 x 12,0V / 5,0Ah 2 x 12,0V / 4,5Ah 2 x 12,0V / 5,0Ah	2 x YB12-5000 2 x MPB12-4500 2 x MPB12-5C

S

STIEGELMEYER

	Batteries compatible with			SYNTHES
	Description / original no.	Туре	Voltage/capacity	axcom product no.
:	Power drive 530.200	NiCd	12,0V / 1,25Ah	MB1370
I	Power battery 532.003	NiCd	12,0V / 0,7Ah	MB1373

IKA	UTW	



Batteries compatible with ...

Description / original no.	Туре	Voltage/capacity	axcom product no.
lift type Elli	lead	12,0V / 12,0Ah	MPB12-12
	lead	12,0V / 10,0Ah	SB51101

Batteries compatible with	Т		TRIMED	
Description / original no.	Ту	rpe Volt	tage/capacity	axcom product no.
MEDICAL AC TESTER (II)	Ni	Cd 10,8	3V / 3,0Ah	MS9200

Batteries compatible with			VÖLKER
Description / original no.	Туре	Voltage/capacity	axcom product no.
2080 / 3080	9 volt e-block	9,0V	DU1604
S960 (NEW version) , 4 batteries	lead	24,0V / 1,2Ah	MB3654
S960-1 (Version A/B/C)	lead	24,0V / 5,0Ah	MB3805
patient bed S960	lead	6,0V / 5,0Ah	MB3655
patient bed S960-1 (Version A/B/C), 2 batteries	lead	24,0V / 7,2Ah	MB3656
patient bed S960-2	lead	12,0/3,4Ah	MB3657

Batteries compatible with			VICTOR
Description / original no.	Туре	Voltage/capacity	axcom product no.
lift Handimove	lead	2 x 12,0 V / 7,2 Ah	2 x MPB12-7200

Batteries	compatible	with



0,

Description / original no.	Туре	Voltage/capacity	axcom product no.
Accuvac suction pump	lead	12,0V / 7,2Ah	GB1270
Accuvac rescue basic	lead	12,0V / 3,4Ah	GB1226
Accuvac rescue	NiCd	12V / 3,0Ah	MB1015
Suction pump accuvac rescue - original	NiCd	12,0V / 3,0Ah	MB1015-O



Batteries compatible with		WISSNER-I	BOSSERHOF
Description / original no.	Туре	Voltage/capacity	axcom product no.
Multicare 2 batteries, je 151 x 65 x 102 mm	lead	24,0V / 7,2Ah	MB3658
Eleganza(Q7)/Estetic 2 batteries, je 97 x 43 x 53 mm	lead	24,0V / 1,2Ah	MB3659

	Batteries compatible with		WE	LCH ALLYN	
-	Description / original no.	Туре	Voltage/capacity	axcom product no.	_
1 A	monitor Atlas 622/623	lead	6,0V / 7,2Ah	GB665	S
Ń	monitor Vital Signs LXI	lead	6,0V / 4,0Ah	MB1501-O	т
	monitor 300er series – original	lead	6,0V / 6Ah	MB4500-O	
P	monitor 53NOP -4200-84	lead	6,0V /4,5Ah	MB5300AX	V

WEINMANN

	Batteries compatible with			ZOLL
	Description / original no.	Туре	Voltage/capacity	axcom product no.
	defibrillator AutoPulse, board type 8700-0752-01 -8700-0752-01	Li-Ion	36,3V / 2,5Ah	MB1490-O
11 2	defibrillator NTP 2 - 1400/1600/1700/M-series/SD2000 SD4410	lead	10,0V / 2,5Ah	MB1400
10	defibrillator M-series CCT/M-series PD4410XL smart	lead	10,8V / 4,5Ah	MB1400XL-O
	defibrillator E-series/R-series BLS/ALS/PLUS original with capacity display	Li-Ion	10,8V / 5,8Ah	MB1440-O

# Cell replacement: Fresh cell culture for your batteries

We can swap used-up cells for new, high-quality industrial cells. The casing and securing elements are retained. You get an absolutely new-quality battery with a 12-month guarantee on the fitted cells.

#### The benefits at a glance:

- **::** 100% performance and up to 60% cost savings compared with original batteries
- :: We only use quality cells and we apply comprehensive quality checks
- :: Certified to ISO 9001:2008
- :: It's easy to do: Send in your original battery and we take care of the rest.
- :: Avoid waste and protect resources



# We offer cell replacements for these batteries:



#### ARJO

Modell.	Туре	voltage / capacity	axcom Product no.
Marisa / Sara3000 / NDA0100-03 / NDA0100-20	lead	24,0V / 4,0Ah	MBU0100
lift Tempo / Opera type KPA0100-0510	lead	24,0V / 5,0Ah	MBU0101
patient lift Trixie / Sarita type KTA0102	lead	24,0V / 2,0Ah	MBU0102
patient lift MAXIMOVE NDA0100-03	lead	24,0V / 4,0Ah	MBU1301
Carendo chair / Maxi Twin lift battery type NAE0100-083	lead	24,0V / 2,5Ah	MBU083
lift type 102193 und Meditechnik bath lift	NimH	12,0V / 7,2Ah	MBB13821 / 849
lift type Carina	lead	24,0V / 2,5Ah	MB9001

### **AKS / ERGOLET**

Modell.		voltage / capacity	axcom Product no.
patient lift FAUN	lead	36,0V / 2,3Ah	MB9001

#### ERGOLET

Modell.		voltage / capacity	axcom Product no.
patient lift Faun	lead	36,0V / 2,3Ah	MB9001
for Ergolet patient lift A/S Modell Maya 791	lead	36,0V / 2,0Ah	MBU791

## HORCHER

Modell.		voltage / capacity	axcom Product no.
patient lift Lexa / Raisa / Diana /Caesar / Ekopack original -No.: L103-1020	lead	24,0V / 2,9Ah	MBU1060
patient lift Nina / Mona original -No.: N206-1012	lead	36,0V / 2,0Ah	MB9001
patient lift Nina / Mona original -No.: N200-1012	lead	36,0V / 2,2Ah	MB9001

#### **MOLIFT**

Modell.		voltage / capacity	axcom Product no.
Molift Power Pac type 1341000	NimH	26,4V / 3,0Ah	MB9001
Molift Power Pac type 0540000	NiCd	14,4V / 1,9Ah	MB9001
Molift Power Pac type 0541000-2	NimH	14,4V / 3,0Ah	MB9001
Molift Power Pac type 1221115	NiCd	26,4V / 3,0Ah	MB9001

## HILL-ROM / LIKO


Modell.		voltage / capacity	axcom Product no.
MIR / Viking M24	NimH	24,0V / 2,5Ah	MB9001
LIKO / Golvo 8000/ 8008 type 2006107 / Linak BAI 20001-02	NimH	24,0V / 2,5Ah	MBU8000

#### (IF SWAPPING CELLS PLEASE SEND IN BATTERY INCLUDING CASING)

# SpO<sub>2</sub>-sensors, ECG-cables and Blood pressure cuff about Axcom.

Sensitive finger and ear sensors, ECG cables: and perfect plug-fit. Make use of Axcom's com-Our products are 100% compatible with and structurally identical to the originals. We guarantee excellent workmanship, high durability

prehensive services, advice and guarantees.



Product no..: MCA-1NE-NE-L SpO2 extension cable for Nellcor 2,4m length

## PHILIPS

Topseller





#### Product no..: MCS-1PH03-L

SpO<sub>2</sub> Soft finger sensor for Philips, such as M1191A/ M1191B Suitable for adults (from 40kg upwards) from 40kg. 3m length. plug-in connection 052/060





#### Product no..: MCS-1PH03-S SpO<sub>2</sub> Soft finger sensor for Philips, such as M1194A Suitable for adults (from 40kg upwards) from 40kg. 90cm length





#### Product no..: MCS-1PH08-L SpO<sub>2</sub> Ohrsensor for Philips, such as M1194A Suitable for adults (from 40kg upwards) from 40kg 3m length. plug-in connection 052/060





#### Product no..: MCS-1PH01-L

SpO<sub>2</sub> finger sensor clip for Philips, such as M1196A (M1191A Softclip). Suitable for adults (from 40kg upwards) from 40kg 3m length





#### Product no..: MCS-2PHS05-L SpO<sub>2</sub> Soft finger sensor for Philips, such as M1192A Suitable for children/small adults. 3m length



#### Product no..: MCA-1PH-PH-L

SpO<sub>2</sub> extension cable for Philips, such as M1941A. suitable for M2601A, M3000A, M3500B 2,4m length

## NIHON KOHDEN



Product no..: MCS-2NIS03-L

3m length



••••

#### MCS-1NI01-S

SpO<sub>2</sub> finger sensor clip for Nihon Kohden Suitable for adults (from 40kg upwards) from 40kg 90cm length



SpO<sub>2</sub> finger sensor Soft for Nihon Kohden, such as TL-201T. Suitable for adults (from 40kg upwards)

**Product no.:: MCA-1NE-NE-L** SpO<sub>2</sub> extension cable for Nellcor 2,4m length





suitable for	Illustration	discription	product no.
Philips 78352A/C, 78354A/C,M1008A/B M1176A,M1350B,M3000A M3001A (with M1598B, M1599B, hoses) Colin , Datascope:	F	<b>Blood pressure cuff Comfort for thighs, such as M1576A,</b> 1 tube for adults with extremities about 42-54cm, latex-free, PVC-free, reusable with replaceable bladder	PCT1T000Y
	G	<b>Blood pressure cuff Comfort, such as M1574A,</b> 1 tube for adults with extremities about 27-35 cm, latex-free, reusable, with replaceable bladder	PCT1A000Y
Passport, Acutor Fukada Denshi	Η	<b>Blood pressure cuff Comfort, such as M1573A,</b> 1 tube for small adults with extremities about 20,5-28cm, latex-free, PVC-free, reusable with replaceable bladder	PCT1S000Y
Spacelabs : all types older Welch-Allyn models with Luer-typee Connector	I	<b>Blood pressure cuff Comfort, such as M1571A,</b> 1 tube for baby 10-15cm, latex-free, PVC-free, reusable with replaceable bladder	PCT11000Y
Criticare, Siemens with Bajonett-connector	A	Blood pressure cuff for thighs, 1 tube for adults with extremities about 42-50cm, latex-free, PVC-free reusable without bladder	PCT1TR000Y
Mindray Goldway	В	<b>Blood pressure cuff,</b> 1 tube for adults with extremities about 32-42cm, latex-free, PVC-free, reusable without bladder	PCT1ALR000Y
	С	<b>Blood pressure cuff,</b> 1 tube for adults with extremities about 24-32cm, latex-free, PVC-free, reusable without bladder	PCT1AR000Y
	D	<b>Blood pressure cuff,</b> 1 tube for adults with extremities about 17-25cm, latex-free, PVC-free, reusable without bladder	PCT1ASR000Y
	E	<b>Blood pressure cuff,</b> 1 tube for children with extremities about 15-22cm, latex-free, PVC-free, reusable without bladder	PCT1CR000Y



## • highest quality

- improved comfort
- reusable
- latex-free
- PVC-free
- replaceable bladder

suitable for	Illustration	discription	product no.
GE-Marquette Datex-Ohmeda	F	<b>Blood pressure cuff Comfort for thighs, such as 5082-88-4</b> 2 tube for adults with extremities about 42-54cm, latex-free, PVC-free, reusable with replaceable bladder	РСТ2Т000Ү
GE(Marquette) with Male quick-connect quick connect, Dinama plus,	G	<b>Blood pressure cuff Comfort, such as 5082-86-4</b> 2 tube for adults with extremities about 27-35 cm latex-free, reusable with replaceable bladder	PCT2A000Y
XL,SX/XP,8100 series, all GE corometrics models	Н	<b>Blood pressure cuff Comfort, such as 5082-85-4</b> 2 tube for small adults with extremities about 20,5-28cm latex-free, PVC-free, with replaceable bladder	PCT25000Y
Critikon Datex-Ohmeda with a Male und a Female guick-	I	<b>Blood pressure cuff Comfort, such as 5082-82-4,</b> 2 tube for baby 10-15cm, latex-free, PVC-free, reusable with replaceable bladder	PCT21000Y
connect quick connect New Welch-Allyn Vitaldatenmodelle with	A	<b>Blood pressure cuff for thighs,</b> 2 tube for adults with extremities about 42-50cm, latex-free, PVC-free, reusable without bladder	PCT2TR000Y
Criticon threaded plug	В	<b>Blood pressure cuff,</b> 2 tube ,for adults with extremities about 32-42cm, latex-free, PVC-free, reusable without bladder	PCT2ALR000Y
	С	<b>Blood pressure cuff,</b> 2 tube for adults with extremities about 24-32cm, latex-free, PVC-free, reusable without bladder	PCT2AR000Y
	D	<b>Blood pressure cuff,</b> 2 tube for adults with extremities about 17-25cm, latex-free, PVC-free, reusable without bladder	PCT1ASR000Y
	E	<b>Blood pressure cuff,</b> 2 tube for children with extremities about 15-22cm, latex-free, PVC-free, reusable without bladder	PCT2CR000Y

# ECG-cables

discription













	ECG master cable for Philips, such as M1669A 3 - core, 2,70m, straight 12-plug	M1001A, M1001B, M1002A, M1002B, M3000A, M3001A, M3002A, M8105A, M8102A, M4735A, M3535A, M3536A, M3536M, 863063, 863064, 78352C, 78354C, 863065, 863066, 863068, 863077, M4735A	M1671A, M1672A (ECG-I3C040GDE), M1673A, M1674A, M1675A, M1678A, M1622A, M1624A, M1678A, M1626A	ECG-IA3370GDE
•	ECG patient lead for Philips, such as M1672A 3 - core with clip - length 102cm. alternativ with press stud	master cable M1669A (ECG-IA- 3370GDE)		ECG-I3C040GDE
	<b>ECG-patient lead</b> 5 times cables with DIN plug, press stud - 61 cm (ro, gr, br, we, sw)	Biosys, Datascope, Nihon Koh- den, AIR SHIELDS	suitable for all DIN-connec- tor-Systeme, Code: C	ECG-EC024S51
	ECG master cable for Philips, such as M1510A 3 - core 2,70m, IEC 12 Pin connector	78352C, 7834C, 78834C (MCS- 2PH01-L), M1001A/B, M1002A/B, M1722A/B (MCS-2PH03-L), M1723A/B (MCS-2PH03-L), M2475B (MCS-2PH03-L), M3000A, M3001A, M4735A, VM4/V24	M1671A, M1672A (ECG-I3C040GDE), M1673A, M1674A, M1675A, M1678A, M1622A, M1624A, M1678A, M1626A	ECG-I3320DE
	ECG patient lead for Philips, such as M1613A 3 - core 1m; IEC clip	M3000A, 862474, 862478, M1001B, M1002B, 863077, 862475, 862483, 863063, 863064, 863065, 863066, 863068, M3002A, M8102A, M8105A, ED040C5A	M1510A (ECG-I3320DE), M1550C	ECG-I3C040DE
	ECG master cable for Philips, such as M1590A 3 - core, 2,7m 8 plug	M1722A, M1722B, M1723A, M1723B, M2475B, M3500B, M5500B, M1649A	M1611A, M1613A (ECG-I3C- 040DE), M1615A, M1619A	ECG-I3-530DE
	ECG patient lead for Philips, such as ED040S3I, 3 - core, 60cm, SHP ACC GRBR 3-core, IEC1 (EUROPA)	M1615A, ED040C5A		ECG-I3S040DE
	ECG patient lead for Dräger, such as EJ024C3I, 3 - core, 60cm, SHP ACC GRBR 3-core, IEC1 (EUROPA)	EJ040S3I, EJ040S5I		ECG-I3C420JE

compatible with

use with cables

Art. No.

	discription	compatible with	use with cables	Art. No.
	ECG master cable for Philips, such as M1530A 5 - core, 2,70m, IEC 12-Pin-connector	M1001A, M1001B, M1002A, M1002B, M3000A, M3001A, M3002A, 863077, 863063, 863064, 863065, 863068, 863066, M3535A, M3536A, M3536M, M4735A, M1722A, M1722B, M1723A, M1723B, M1724A, M2475B, M1649A, M8102A, M8105A, 78352C, 7834C, 78834C, M3000A, M3001A, VM4/V24	CMS-electrode cable; 5 - core	ECG-15320DE
	ECG patient lead for Philips, such as M1633A 5 - core, 1,0m, IEC clip, LA,RA,V	master cable M1530A (ECG- I5320DE)		ECG-15C040DE
	ECG master cable for Philips, such as M1663A 10 - core, 2,7m	M1663A, M1949A		ECG-NC040GDV
	ECG patient lead for Philips, such as M1971A 5 - core, 1,0m, IEC-clip, LA,RA,V	master cable M1663A (ECG- NC040GDV); M1949A; M1978A (ECG-INC040GDV)		ECG-15C040GDE
(I)	ECG patient lead for Dräger, such as EJ040C5I, 5 - core, 1,0m, IEC1, 5 - core lead wire, standard	EJ024C3I (ECG-I3C420JE)		ECG-I5C040JE
	ECG master cable for Philips, such as M1668A 5 - core, 2,7m	M1001A, M1001B, M1002A, M1002B, M3000A, M3001A, M3002A, M8105A, M8102A, M4735A, M3535A, M3536A, M3536M, 78352C, 78354C, 860335, 863064, 863065, 863066, 863068, M4735A, M1649A	M1968A, M1644A, M1973A, M1647A, M1971A (ECG-I5C040GDE), M1645A, M1974A, M1648A	ECG-IA5370GDE
and a	ECG-ECG multimed cables for Dräger/Siemens, such as 3368391 / 5950196, 2,4m, with disinfection-resistant color markings	SC6000, SC6002XL, SC7000, SC9000	M1603A	ECG-IA5650QE
	ECG master cable for Hellige/ GE/Marquette, such as 223 418 08, 10 polig 2,2m, IEC, 15 Pin connector	MAC500/1100/1200		ECG-IAB800EV
	electrode lead for Cardiosmart, such as GE 38401817, set of 10 4mm banana plug	Cardiosmart & MAC		ECG-INB800EV
	<b>ECG patient lead for Philips,</b> <b>such as M1978A</b> 2 x 5 core, chest, clip 1 m for master cable M1668A and M1663A	M1668A, M1663A		ECG-INC040GDV

# Certified performance for every medical equipment application

Axcom batteries for medical equipment comply with the stringent demands of the MPG (Medizinproduktgesetz, Germany's Medical Devices Act). The company's own safety officer makes sure of this, pursuant to ISO 9001:2008.

#### 360° reliability

Years of experience show you can rely on Axcom batteries when it matters most.

- :: Injection and infusion pumps :: Patient monitors
- :: Defibrillators

:: Patient beds

:: Patient lifts

### Top quality – and that includes our service & support

For many years Axcom has been a trusted partner in the medical equipment sector, offering a comprehensive range of high-end batteries and accessories as well as professional support by specially trained consultants. Our engineers and technicians are available on the phone at any time and will gladly come to your premises to present you with the best possible solutions for your application area.

# More power - more efficiency: Primary batteries and individual cells.

#### Primary batteries and individual cells

Our range also includes non-rechargeable primary batteries including ones made by Duracell, Varta, Energizer and Camelion, for use in devices such as pagers. We also offer rechargeable individual cells by eneloop, Panasonic, Sanyo and other big-name manufacturers.



Product no.: Z-YT3000SCP NiMh battery SUB C/YT3000SCP 1.2V, 3Ah



Product no.: Z-HR-3UTGB NiMh battery PANASONIC eneloop LR06/MIGNON/AA, 1.2V, 2Ah



Product no.: Z-HR-3U-2700 NiMh battery **Fujitsu** (such as Sanyo HR3U) HR-3U LR06 MIGNON/AA, 1.2V 2.7Ah

# Custom solutions manufactured by us.

Do you need specialized batteries or battery packs for your application? Please tell us your specifications or let us determine them by talking with you. Our engineers can configure the battery which suits your needs precisely. You will be given detailed data sheets and, if necessary, approval samples before production begins. You can be sure that you have the right partner at your side in terms of quality, punctuality and value for money.





## Visit our newest shop on the web.



Tidy, clear and up-to-date: In our completely redesigned online shop, you can quickly find your desired product - it's worth the visit.

#### **Our products**

Batteries for medical equipment, accessories for medical equipment, batteries for radio devices and pagers, accessories for devices, camera batteries, lamp batteries, controller batteries, batteries for electrical tools, scanner batteries, scanner accessories, batteries for mobile printers, radar measurement device batteries, lead batteries, charging equipment, individual cells, hearing aid batteries, primary batteries, button cells, customized products, cell swapping

#### Our service

- :: Structured organization of product management, procurement, technology, quality control, marketing, sales and logistics with the constant aim of providing our customers with the best products at the right time.
- :: Our committed sales team is always up to date because it is so widely represented in the market, and it gladly shares its knowledge with our customers.
- :: Our engineers and technicians know the latest developments in battery technology. This means our products are always a step ahead technically.
- :: Our technicians maintain direct contact with our customers by phone and on site, which means they know their specific requirements and can offer the best possible solutions.
- :: Our marketing activities and sales documentation help our customers with their own marketing.

#### Offer service

:: We can draw up an offer specially for you which takes your wishes into account.

#### **Efficient delivery**

:: Because we have so many products in stock, we can deliver most of our products anywhere in Germany within 24 hours. We will gladly tell you how long it will take to deliver products to the rest of Europe or anywhere else in the world.

#### Axcom can provide you with the ultimate quality accessories for your equipment.

#### Contact

Tel.: +49 (0) 21 54 - 48 38-0 or 0800 - 253 98 35 Fax: +49 (0) 21 54 - 48 38-444 or 0800 - 253 93 29



11.2018



Our General Terms and Conditions of Business apply. All of our products are alternative products and replicas we have built ourselves, and are not the original products made by the manufacturers of devices. Some of the brand and company names used here are registered trademarks and are used solely to indicate compatibility. We shall not be liable for mistakes or printing errors.

