


# MFM-CNS

Central Monitoring System



**MFM-CNS**  
Central Monitoring System



## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



A world of potential

Edan Instruments, Inc. | 3/F-B, Nanshan Medical Equipments Park,  
1019# Nanhai Rd., Shenzhen | 518067 P.R. China  
+86.755.26898326 | [www.edan.com.cn](http://www.edan.com.cn) | [info@edan.com.cn](mailto:info@edan.com.cn)



ENG-OB & GYN-MFM-CNS-V3.0-20150112

© 2015 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.



A world of potential

## Web Viewer

The web viewer function helps doctors with remote monitoring from their own PC, improving communication efficiency.



F9 Express



F9



F6



F3



F2



Single or dual screen displays



Wi-Fi / LAN Connctions



OB Office



Nurse Station



Labor & Delivery Room



Login by appointed user authorization

Username

CNS is designed different extent of authority of users

Alerting status indicator

Alert sourcing indicator



Device indicators


Department allocations

Fetal monitoring

Fetal and maternal monitoring


Maternal monitoring

### CTG Analysis



To assist diagnosing , diagnostic result can be saved in electronic format.

### Partogram



To provide an accurate record of the progress in labor, so that any delays or deviations from normals can be presented immediately and treated accordingly.

# MFM-CNS

## Central Monitoring System

MFM-CNS is a comprehensive central nursing system for the obstetric cares. It is an integrated system that covers the entire continuum of obstetric care across the whole pregnancy stage from the first antepartum visit to post delivery care. The system displays data and information transmitted from fetal and maternal monitors and generates alerts for critical events. Easy access to well-organized information improves efficiency, allowing you to focus more attention on patients and achieve the best quality of care.

- Web viewer for remote monitoring from doctors' PCs
- Built-in antepartum CTG analysis
- Display up to 6 waveforms of each bedside monitor simultaneously
- Review of events, CTG waveforms, trends and NIBP data
- Large capacity of storage up to 100,000 records
- Wired & wireless central monitoring system solution
- Connection of up to 32 bedside fetal/maternal monitors
- Up to 24 hours monitoring data for each record
- Up to 200 groups of NIBP data
- Export PDF format report