



F9 Series
Fetal & Maternal Monitor

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | 3/F-B, Nanshan Medical Equipments Park,
1019# Nanhai Rd., Shenzhen | 518067 P.R. China
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn

© 2015 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.



ENG-OB & GYN-F9-V3.0-20150112

F9 Series

Fetal & Maternal Monitor



F9 Series


Fetal & Maternal Monitor

Quick setup and simplified workflow

The 'START' button could be configured to integrate patient info and printing. This could help the operator to simplify the workflow and work with one button only.



Built-in antepartum CTG analysis



The antepartum CTG analysis provides objective assessment of CTG for the caregivers' reference and ensures easy access to read FHR patterns.

Various printing options

F9 series are compatible with 152 and 150 mm thermal printing paper in either American or International standard.

Telemetry transducer system

As an optional accessory of F9, the telemetry transducer system, with its lightweight and waterproof probes offers wireless monitoring with great freedom and convenience, especially for water delivery.

60 hours memory for seamless monitoring



12.1" color touch screen

The touch-screen, color display allows quick setup and optimized performance.

F9 series fetal & maternal monitor provide solutions for all areas of fetal monitoring, ranging from antepartum, intrapartum, and postpartum clinical applications. With integrated monitoring of twins' FHR, uterine activity, fetal movement, intrauterine pressure and direct ECG as well as maternal NIBP, SpO₂, ECG and TEMP, F9 series offer the most cost effective and flexible solution for all your fetal monitoring needs, without compromising quality, precision, performance and ergonomics.

Multi-display Modes



Switching among three display modes, user may choose the most suitable one for clinical use.

12-crystal and 1MHz waterproof transducer

| Depth | 7-crystal | 9-crystal | 12-crystal |
|-------|-----------|-----------|------------|
| 40mm | | | |
| 80mm | | | |
| 120mm | | | |

Our 12-crystal ultrasound probe employs advanced technology unique to EDAN, providing wider beam area and more homogeneous signal, thus realizing better performance for bedside monitors than conventional transducers.