



Catalogue Optiques rigides

Immed Europe
Parc d'activité de le Broye
59710 Ennevelin
Tél : 03.28.55.02.55 Fax : 03.28.55.02.58

Immed Benelux
Avenue Winston Churchill 192A
1180 Bruxelles
Tél : 02.381.15.55 – Fax : 02.381.15.61



Arthroscope

Tous nos arthroscopes sont autoclavables 18min à 134°C



Art n°	Désignation
10-0001-00	Arthroscope HD Diam :2.7mm - long :110mm – ang : 0°
10-0002-00	Arthroscope HD Diam :2.7mm - long :110mm – ang : 30
10-0003-00	Arthroscope HD Diam :2.7mm - long :110mm – ang : 70°
10-0004-00	Arthroscope HD Diam :2.7mm - long :175mm – ang : 0°
10-0005-00	Arthroscope HD Diam :2.7mm - long :175mm – ang : 30°
10-0368-00	Arthroscope HD Diam :2.7mm - long :175mm – ang : 45°
10-0006-00	Arthroscope HD Diam :2.7mm - long :175mm – ang : 70°
10-0007-00	Arthroscope HD Diam :2.7mm - long :187.5mm – ang : 0°
10-0008-00	Arthroscope HD Diam :2.7mm - long :187.5mm – ang : 30°
10-0009-00	Arthroscope HD Diam :2.7mm - long :187.5mm – ang : 70°
10-0163-00	Arthroscope HD Diam :2.7mm - long :141mm – ang : 0°
10-0164-00	Arthroscope HD Diam :2.7mm - long :141mm – ang : 30°
10-0414-00	Arthroscope HD Diam :2.7mm - long :141mm – ang : 45°
10-0165-00	Arthroscope HD Diam :2.7mm - long :141mm – ang : 70°
10-0010-00	Arthroscope HD Diam:4mm - long: 175mm - ang: 0°
10-0011-00	Arthroscope HD Diam :4mm - long :175mm - ang : 30°
10-0085-00	Arthroscope HD Diam :4mm - long :175mm - ang : 45°
10-0012-00	Arthroscope HD Diam :4mm - long :175mm - ang : 70°
10-0223-00	Arthroscope HD Diam :3mm - long :110mm - ang : 0°
10-0224-00	Arthroscope HD Diam :3mm - long :110mm - ang : 30°
10-0261-00	Arthroscope HD Diam :2.4mm - long :98mm - ang : 0°
10-0262-00	Arthroscope HD Diam :2.4mm - long :98mm- ang : 30°
10-0339-00	Arthroscope HD Diam :1.9mm - long :69mm - ang : 0°
10-0340-00	Arthroscope HD Diam :1.9mm - long :69mm - ang : 30°
10-0341-00	Arthroscope HD Diam :1.9mm - long :102mm - ang : 0°
10-0342-00	Arthroscope HD Diam :1.9mm - long :102mm - ang : 30°

Otoscope

Autoclavables 18min à 134°C



Art n°	Désignation
10-0032-00	Otoscope HD Diam :2.7mm – long :50mm – ang : 0°
10-0033-00	Otoscope HD Diam :2.7mm – long :50mm – ang : 30°
10-0034-00	Otoscope HD Diam :4mm – long :50mm – ang : 0°
10-0035-00	Otoscope HD Diam :4mm – long :50mm – ang : 30°
	Otoscope avec speculum d'oreille
10-0134-00	Otoscope HD Diam :4mm – long :50mm – ang : 0°



Sinuscope

Tous nos sinuscopes sont autoclavables 18min à 134°C



Art n°	Désignation
10-0036-00	Sinuscope HD Diam :2.7mm - long :110mm - ang : 0°
10-0037-00	Sinuscope HD Diam :2.7mm - long :110mm - ang : 30°
10-0418-00	Sinuscope HD Diam :2.7mm - long :110mm - ang : 45°
10-0038-00	Sinuscope HD Diam :2.7mm - long :110mm - ang : 70°
10-0039-00	Sinuscope HD Diam :2.7mm - long :175mm - ang : 0°
10-0040-00	Sinuscope HD Diam :2.7mm - long :175mm - ang : 30°
10-0387-00	Sinuscope HD Diam :2.7mm - long :175mm - ang : 45°
10-0041-00	Sinuscope HD Diam :2.7mm - long :175mm - ang : 70°
10-0042-00	Sinuscope HD Diam :2.7mm - long :187.5mm - ang : 0°
10-0043-00	Sinuscope HD Diam :2.7mm - long :187.5mm - ang : 30°
10-0044-00	Sinuscope HD Diam :2.7mm - long :187.5mm - ang : 70°
10-0466-00	Sinuscope HD Diam :2.7mm - long :141mm - ang : 0°
10-0467-00	Sinuscope HD Diam :2.7mm - long :141mm - ang : 30°
10-0486-00	Sinuscope HD Diam :2.7mm - long :141mm - ang : 45°
10-0045-00	Sinuscope HD Diam :4mm - long :175mm - ang : 0°
10-0046-00	Sinuscope HD Diam :4mm - long :175mm - ang : 30°
10-0086-00	Sinuscope HD Diam :4mm - long :175mm - ang : 45°
10-0047-00	Sinuscope HD Diam :4mm - long :175mm - ang : 70°
10-0474-00	Sinuscope HD Diam :4mm - long :175mm - ang : 90°
10-0264-00	Sinuscope HD Diam :2.4mm - long :98mm - ang : 0°
10-0265-00	Sinuscope HD Diam :2.4mm - long :98mm - ang : 30°
10-0416-00	Sinuscope HD Diam :1.9mm - long :69mm - ang : 30°
10-0422-00	Sinuscope HD Diam :1.9mm - long :69mm - ang : 0°
10-0573-00	Sinuscope HD Diam :3mm - long :110mm - ang : 0°
10-0574-00	Sinuscope HD Diam :3mm - long :110mm - ang : 30°

Sinuscope Grand Angle 180°



10-0426-00	Sinuscope HD 180° Diam :4mm – long:175mm – ang: 30°
10-0622-00	Sinuscope HD 180° Diam :4mm – long:175mm – ang: 45°

Microlaryngoscope



Autoclavables 18min à 134°C

Art n°	Désignation
10-0090-00	Micro Laryngoscope HD Diam :4mm – long:175mm – ang: 70°
10-0219-00	Micro Laryngoscope HD Diam :4mm – long:175mm – ang: 90°
10-0475-00	Micro Laryngoscope HD Diam :5.8mm – long:188mm – ang:70°




Laryngoscope

Autoclavables 18min 134°C

Art n°	Désignation
10-0177-00	Laryngoscope HD Diam :8mm – long :190mm – ang: 70°
Laryngoscope "illumination   latérale "	
10-0476-00	Laryngoscope HD Diam :10mm – long:190mm – ang: 70°
10-0477-00	Larungoscope HD Diam :10mm – long:190mm – ang: 90°

Telelaryngoscope

Nos telelaryngoscopes sont immergeables



Art n°	Désignation
10-0042-00	Tele Laryngoscope HD Diam :11mm – long:200mm – ang: 90°
10-0468-00	Tele Laryngoscope HD Diam :11mm – long:200mm – ang: 70°

Bronchoscope

Autoclavables 18min à 134°C

Art n°	Désignation
10-0050-00	Bronchoscope HD Diam:5.5mm – long:490mm – ang: 0°
10-0051-00	Bronchoscope HD Diam:5.5mm – long:490mm – ang: 30°
10-0383-00	Bronchoscope HD Diam:2.9mm – long:365mm – ang: 0°
10-0384-00	Bronchoscope HD Diam:2.9mm – long:365mm – ang: 30°
10-0618-00	Bronchoscope HD Diam:2.9mm – long:303mm – ang: 0°
10-0619-00	Bronchoscope HD Diam:2.9mm – long:303mm – ang: 30°

Oesophagoscope

Autoclavables 18min à 134°C

Art n°	Désignation
10-0049-00	Oesophagoscope HD Diam:5.5mm - long:530mm – ang: 0°




Hystéroscope

Autoclavables 18min à 134°C

Art n°	Désignation
10-0018-00	Hystéroscope HD Diam:2.7mm – long:302mm – ang: 0°
10-0019-00	Hystéroscope HD Diam:2.7mm – long:302mm – ang: 12°
10-0162-00	Hystéroscope HD Diam:2.7mm – long:302mm – ang: 30°
10-0020-00	Hystéroscope HD Diam:2.7mm – long:302mm – ang: 70°
10-0021-00	Hystéroscope HD Diam:4mm – long:302mm – ang: 0°
10-0022-00	Hystéroscope HD Diam:4mm – long:302mm – ang: 12°
10-0023-00	Hystéroscope HD Diam:4mm – long:302mm – ang: 30°
10-0024-00	Hystéroscope HD Diam:4mm – long:302mm – ang: 70°
10-0395-00	Hystéroscope HD Diam:2.9mm – long:302mm – ang: 0°
10-0393-00	Hystéroscope HD Diam:2.9mm – long:302mm – ang: 12°
10-0385-00	Hystéroscope HD Diam:2.9mm – long:302mm – ang: 30°
10-0615-00	Hystéroscope HD Diam:2mm – long:260mm – ang: 0°
10-0617-00	Hystéroscope HD Diam:2mm – long:260mm – ang: 12°
10-0616-00	Hystéroscope HD Diam:2mm – long:260mm – ang: 30°

Laparoscope 90/90 – coudé 2x90


Autoclavable 18min à 134°C



Art n°	Désignation
10-0138-00	Laparoscope HD 90/90 – coudé 2x90 Diam:10mm – long:268mm – ang : 0°, 1 canal de travail – Diam:6mm – long:329mm En option: avec adaptateur réf n° 24-0092-00 – canal de travail long:372mm

Laparoscope

Autoclavable 18min à 134°C



Art n°	Désignation
10-0013-00	Laparoscope HD Diam:5mm – long:300mm – ang: 0°
10-0014-00	Laparoscope HD Diam:5mm – long:300mm – ang: 30°
10-0109-00	Laparoscope HD Diam:5mm – long:300mm – ang: 45°
10-0225-00	Laparoscope HD Diam:5mm – long:360mm – ang: 0°
10-0226-00	Laparoscope HD Diam:5mm – long:360mm – ang: 30°
10-0227-00	Laparoscope HD Diam:5mm – long:360mm – ang: 45°
10-0326-00	Laparoscope HD Diam:3mm – long:300mm – ang: 30°
10-0327-00	Laparoscope HD Diam:3mm – long:300mm – ang: 0°
10-0390-00	Laparoscope HD Diam:5mm – long:450mm – ang: 0°
10-0472-00	Laparoscope HD Diam:5mm – long:450mm – ang: 30°



Laparoscope “plein écran”

Art n°	Désignation
10-0483-00	Laparoscope HD Diam:5.5mm – long:310mm – ang: 0°
10-0484-00	Laparoscope HD Diam:5.5mm – long:310mm – ang: 30°
10-0063-00	Laparoscope HD Diam:6.5mm – long:310mm – ang: 0°
10-0064-00	Laparoscope HD Diam:6.5mm – long:310mm – ang: 30°
10-0596-00	Laparoscope HD Diam:6.5mm – long:350mm – ang: 0°
10-0597-00	Laparoscope HD Diam:6.5mm – long:350mm – ang: 30°


Laparoscope – avec filtre de chaleur

Autoclavables 18min à 134°C

Art n°	Désignation
10-0317-00	Laparoscope HD Diam:10mm – long:310mm – ang: 45°
10-0329-00	Laparoscope HD Diam:10mm – long:450mm – ang: 0°
10-0330-00	Laparoscope HD Diam:10mm – long:450mm – ang: 30°
10-0331-00	Laparoscope HD Diam:10mm – long:450mm – ang: 45°
10-2315-00	Laparoscope HD Diam:10mm – long:310mm – ang: 0°
10-2316-00	Laparoscope HD Diam:10mm – long:310mm – ang: 30°

Laparoscope “plein écran”

Autoclavable 18min à 134°C




Art n°	Désignation
10-0089-00	Laparoscope HD Diam:10mm – long:310mm – ang: 45°
10-0229-00	Laparoscope HD Diam:10mm – long:450mm – ang: 0°
10-0230-00	Laparoscope HD Diam:10mm – long:450mm – ang: 30°
10-0231-00	Laparoscope HD Diam:10mm – long:450mm – ang: 45°
10-0232-00	Laparoscope HD Diam:10mm – long:540mm – ang: 0°
10-0233-00	Laparoscope HD Diam:10mm – long:540mm – ang: 30°
10-2015-00	Laparoscope HD Diam:10mm – long:310mm – ang: 0°
10-2016-00	Laparoscope HD Diam:10mm – long:310mm – ang: 30°



Thoracoscope

Autoclavable 18min à 134°C



Art n°	Désignation
10-0065-00	Thoracoscope HD Diam:10mm – long:210mm – ang: 0°
10-0066-00	Thoracoscope HD Diam:10mm – long:210mm – ang: 30°

Cystoscope

Autoclavable 18min à 134°C



Art n°	Désignation
10-0025-00	Cystoscope HD Diam:2.7mm – long:302mm – ang: 0°
10-0139-00	Cystoscope HD Diam:2.7mm – long:302mm – ang: 12°
10-0026-00	Cystoscope HD Diam:2.7mm – long:302mm – ang: 30°
10-0027-00	Cystoscope HD Diam:2.7mm – long:302mm – ang: 70°
10-0028-00	Cystoscope HD Diam:4mm – long:302mm – ang: 0°
10-0029-00	Cystoscope HD Diam:4mm – long:302mm – ang: 12°
10-0030-00	Cystoscope HD Diam:4mm – long:302mm – ang: 30°
10-0031-00	Cystoscope HD Diam:4mm – long:302mm – ang: 70°
10-0343-00	Cystoscope HD Diam:1.9mm – long:187.5mm – ang: 30°
10-0359-00	Cystoscope HD Diam:1.9mm – long:187.5mm – ang: 0°
10-0347-00	Cystoscope HD Diam:2.7mm – long:187.5mm – ang: 0°
10-0348-00	Cystoscope HD Diam:2.7mm – long:187.5mm – ang: 30°
10-0351-00	Cystoscope HD Diam:2.7mm – long:187.5mm – ang: 70°
10-0386-00	Cystoscope HD Diam:2.9mm – long:302mm – ang: 30°
10-0394-00	Cystoscope HD Diam:2.9mm – long:302mm – ang: 12°
10-0396-00	Cystoscope HD Diam:2.9mm – long:302mm – ang: 0°
10-0632-00	Cystoscope HD Diam:2mm – long:260mm – ang: 0°
10-0633-00	Cystoscope HD Diam:2mm – long:260mm – ang: 12°
10-0634-00	Cystoscope HD Diam:2mm – long:260mm – ang: 30°