

# Superior Technology Incredible Performance

## F15

Fetal & Maternal Monitor



F15

F15 Air



### About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Global Headquarters:**  
Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,  
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122 P.R. China  
+86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com.cn](mailto:info@edan.com.cn)

**U.S. and Canada inquiries:**  
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnostics.com](http://www.edandiagnostics.com) | [edan-info@edandiagnostics.com](mailto:edan-info@edandiagnostics.com)

©Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

ENG-OB&GYN-F15-V1.0-20191111  
For Medical use only



# Evolutionary Generation of Fetal Monitoring

As an ultimate CTG machine, F15 series fetal and maternal monitor is the perfect combination of advanced technology and sophisticated algorithm for Obstetrics' new era. With state-of-art design, delicate appearance and strong mechanical structure, F15 presents the midwives with the deepest sincerity to bring the most efficient work flow to clinics, such as automatic probe activation and wireless charging, to keep the monitoring process in hand all the time. Consider F15 not only for today's monitoring needs but also as an enduring platform for the future. There are two models for different clinical environment, the wired and the wireless configurations.

-   
 15 mins Backup Memory
-   
 Monitoring Pause
-   
 Transducer Plug and Play
-   
 Lightweight design
-   
 Triplets Monitoring
-   
 FHR&MHR




# Free for movement

Bring freedom in OB clinics by offering wireless accessories including fetal monitoring and maternal monitoring, midwives are no longer being fussy by messy wires surrounding the equipments and the patients. Even more, F15 supports non-stop monitoring when transferring mom-to-be from labor room to delivery room. By 'TRANSFER' function, patients data could be switch from one device to another one at one click.




# Sophisticated alarm function

The essential purpose of monitoring fetal heart rate is to be alert to abnormal fetal heart rate condition. Since CTG is non-invasive long term monitoring methods, it is difficult to be alert to severe condition timely. To include more sophisticated alarm information on F15 gives more power to midwives for quick identifying critical clinical case.



Sinusoidal fetal heart rate



Low variability fetal heart rate